

MEDIC ASSAULT KIT

Learn more about this and other Products with a Mission® at www.NARescue.com

Middle hinge design allows contouring to the body



A Contoured Fit For Room Breaching And Clearing.

The Medic Assault Rescue Kit is our newest compact, low-profile trauma & rescue aid bag. This bag is unique because it is designed with a middle hinge, allowing a contoured fit for room breaching and clearing.

Two main bag compartments are accessible from each end of bag independently. The front of the bag has two large front pockets with full MOLLE on the exterior. Inside, there are four large mesh inner pouches that are secured to the bag interior with hook & loop attachment and provide organization via removable hook & loop labels for customization. The Basic Rescue kit comes with two Bleed pouches, one Splinting pouch, and one Airway pouch. The Advanced Rescue kit comes with two Bleed pouches, one IV pouch, and one Airway pouch.

MOLLE on back and bottom of the exterior of the nylon bag provides additional carry capability. All silent zipper pulls. Removable backpack straps allow attachment directly to body armor if desired (requires additional straps/hardware – not included).

Available in six 500D Cordura colors, including MTC Alpine.



MEDIC ASSAULT KIT

ITEM #	COLOR	OPTION	NSN
80-1650	Black	Basic	----
80-1651	Coyote	Basic	----
80-1653	OD Green	Basic	----
80-1654	Ranger Green	Basic	----
80-1652	Multicam	Basic	----
80-1649	MTC Alpine	Basic	----
80-1731*	Black	Advanced	----
80-1732*	Coyote	Advanced	----
80-1734*	OD Green	Advanced	----
80-1735*	Ranger Green	Advanced	----
80-1733*	Multicam	Advanced	----
80-1730*	MTC Alpine	Advanced	----

* Advanced kits are **R_X** ONLY

FEATURES/SPECIFICATIONS

- Slim Design and Low Empty Weight
- Middle hinge design to allow contouring to the body decreasing profile
- Large mesh internal bags
- Attachable to body armor
- Available in six 500D Cordura colors, including MTC Alpine

DIMENSIONS





















- L 22 in. x W 11.5 in. x D 5 in. (3 in. seam to seam)
- Advanced Weight: 12.5 lbs



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

MEDIC ASSAULT KIT

all KITS INCLUDE:

 QTY 2 C-A-T® Tourniquet	 QTY 4 NAR® Gauze, Wound Packing	 QTY 2 ETD® 6 in. Dressing	 QTY 2 ETD® 4 in. Dressing	 QTY 1 ETD® Abd./Stump	 QTY 2 NAR® Naso w/ Lube, 28F
 QTY 1 Tactical Suction Device	 QTY 2 HyFin® Vent Chest Seal, Twin Pk	 QTY 2 Emergency Survival Wrap	 QTY 1 Ex-Tract Strap	 QTY 2 SAM® Splint II	 QTY 2 Triangular Bandage
 QTY 2 PES Eye Shield	 QTY 1 NAR® Trauma Shears, Lg	 QTY 1 Scissor Leash™	 QTY 1 Aviation Panel, 35" x 29"	 QTY 2 Combat Casualty Card	
 QTY 1 Nametape AIRWAY	 QTY 2 Nametape BLEED	 QTY 1 Permanent Marker, Large	+plus, by kit type:		 QTY 1 Nametape SPLINTING

KIT CONTENTS

- 2 x C-A-T® Tourniquet - Black
 - 4 x NAR® Wound Packing Gauze - 3 in. x 5 yd.
 - 2 x NAR® Emergency Trauma Dressing (ETD®) - 6 in.
 - 2 x NAR® Flat Emergency Trauma Dressing (ETD®) - 4 in.
 - 1 x Abdominal/Stump Emergency Trauma Dressing (ETD®) - 12 in. x 12 in.
 - 2 x NAR® Nasopharyngeal Airway - 28FR - Pre-lubricated
 - 1 x Tactical Suction Device
 - 2 x HyFin® Vent Chest Seal Dressing (Twin Pack)
 - 2 x NAR® Emergency Survival Wrap
 - 1 x Ex-Tract Strap
 - 2 x SAM® Splint II
 - 2 x Triangular Bandage - OD Green
 - 2 x Polycarbonate Eye Shield (PES) - Orange
 - 1 x NAR® Trauma Shears, Black - Large - 7.25 in.
 - 1 x Trauma Shear Scissor Leash™ - Black
 - 1 x Aviation Panel Kit - 35 in. x 29 in.
 - 2 x Combat Casualty Card
 - 2 x BLEED Nametape - Black/Gold
 - 1 x AIRWAY Nametape - Black/Gold
 - 1 x Permanent Marker, Black - Large
- +plus, by kit type:**
- Basic**

 - 1 x SPLINTING Nametape - Black/Gold

Advanced

 - 2 x Quikclot® Z-Fold Combat Gauze - 3 in. x 4 yd.
 - 1 x Cyclone® Pocket Bag Valve Mask (BVM)
 - 1 x Bougie Aided Crick (BAC) Pack (Cricotrydotomy Kit)
 - 2 x Enhanced ARS™ Decompression Needle - 10 ga x 3.25 in.
 - 2 x Saline Lock Kit
 - 1 x BOA® IV Constricting Band
 - 1 x Intraosseous EZ-IO Manual Needle (Box of 6)
 - 1 x NAR® Mini Duct Tape - 2 in.
 - 1 x Halstead Straight Mosquito Forceps - 5 in.
 - 1 x Sharps Shuttle
 - 1 x Armadillo Medicine Case - Black
 - 1 x IV Nametape - Black/Gold

 QTY 2 Combat Gauze® Z-Fold	 QTY 1 Cyclone® Pocket BVM	 QTY 1 BAC-Pack Bougie-Aided Cric	 QTY 2 Enhanced ARS Needle, 10G	 QTY 2 NAR® Saline Lock Kit	 QTY 1 BOA® IV Constricting Band
 QTY 1 Intraosseous EZ-IO Manual Needle	 QTY 1 NAR® Mini Duct Tape, 2 in.	 QTY 1 5" Halstead Straight Mosquito Forceps	 QTY 1 Sharps Shuttle	 QTY 1 Armadillo Med Case	 QTY 1 Nametape IV

