

OPERATOR BLS/IFAK

The Operator BLS/IFAK™ is designed to equip the individual operator with the Basic Life Support tools needed for immediate Self-Aid/Buddy-Aid. This Improved First Aid Kit is configured to address leading causes of preventable death due to penetrating injuries. The versatile BLS IFAK can be worn by the operator in one of three methods, drop down leg rig, attaché to the assault vest/body armor with MOLLE-style connectors or attached directly to the duty belt.

Learn more about this and other Products with a Mission® at www.NARescue.com



FEATURES/SPECIFICATIONS

- Sand & water resistant, self-repairing zipper with silent 550 cord pulls
- MOLLE Style Attachment System Compatible
- Removable Belt and Leg Straps
- Leg Strap with Quick-Grip™

OPERATOR BLS/IFAK

ITEM #	COLOR	NSN
80-0117	Black	----
80-0118	OD Green	----

* Explore more IFAK's at www.NARescue.com

DIMENSIONS


- Case Closed: H 8.5 in. x W 8.75 in. x D 4.5 in.
- Weight: 1 lb 12 oz



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

OPERATOR BLS/IFAK

all KITS INCLUDE:



KIT CONTENTS

- 1x Bag (IFAK)
- 6x Black Talon® Nitrile Trauma Gloves (3 pr.)
- 1x NAR® Nasopharyngeal Airway with Lube - 28F(117 mm)
- 1x HyFin® Vent Chest Seal Twin Pack
- 1x C-A-T® (Combat Application Tourniquet®)
- 1x ETD® 6 in. Emergency Trauma Dressing
- 2x S-Rolled Gauze™ (4.5 in. x 4.1 yd)
- 1x Polycarbonate Eye Shield
- 1x Permanent Marker
- 1x Combat Casualty Card



OPERATOR BLS/IFAK

ITEM #	COLOR	NSN
80-0117	Black	----
80-0118	OD Green	----

* Explore more IFAK's at www.NARescue.com

Learn more about this and other Products with a Mission® at www.NARescue.com

