

S.T.O.R.M. OPERATOR IFAK

The S.T.O.R.M. Operator IFAK is the official kit for the Specialized Tactics for Operational Rescue and Medicine (S.T.O.R.M.) courses developed through a joint effort between the National Tactical Officers Association and the Medical College of Georgia-Center of Operational Medicine.

The S.T.O.R.M. Operator IFAK is one of our smallest Improved First Aid Kits (IFAK) that was specifically designed to mitigate the three leading causes of preventable death in the tactical environment for penetrating trauma (60% extremity bleeding, 33% tension pneumothorax, 6% airway obstruction). This kit provides the necessary equipment for the user/operator to administer immediate self-aid/buddy-aid to injured casualties at the point of wounding. Utilizing Tactical Combat Casualty Care (TCCC) Guidelines and lessons learned from the battlefield, this kit provides law enforcement officers and other public safety personnel with the needed tools to reduce the unnecessary loss of life due to the preventable causes of death encountered in the tactical environment. Made to meet the demands of the austere environment, our modular system can be worn as a drop down leg rig, attached directly to the duty belt or attached to the assault vest or body armor with MOLLE/PALS-style connectors.



Learn more about this and other Products with a Mission® at www.NARescue.com



FEATURES/SPECIFICATIONS

- SEL#: 09ME-01-BAGM
- One of NAR's most compact IFAK profiles
- Made of durable CORDURA® Nylon
- Sand & water resistant, self-repairing zipper with silent 550 cord pulls
- MOLLE/PALS-style Attachment System Compatible
- Removable Belt and Leg Straps
- Leg Strap with Quick-Grip™

Rx ONLY

S.T.O.R.M. OPERATOR IFAK

ITEM #	COLOR	NSN
80-0133	Black	----
80-0135	Coyote	----
80-0134	OD Green	----

* Explore other IFAK's at www.NARescue.com

DIMENSIONS

- Case Closed: H 7.4 in. x W 6.9 in. x D 5 in.
- Weight: 1 lb 12 oz

S.T.O.R.M. OPERATOR IFAK

all KITS INCLUDE:



KIT CONTENTS

- 1x Bag (STORM)
- 2 pr Black Talon® Nitrile Trauma Gloves
- 1x NAR® Nasopharyngeal Airway with Lube - 28F(117 mm)
- 1x HyFin® Vent Chest Seal Twin Pack
- 1x ARS® Needle Decompression Kit (14 G x 3.25 in.)
- 1x C-A-T® - Black (Combat Application Tourniquet®)
- 1x ETD 6 in. Emergency Trauma Dressing
- 2 x S-Rolled Gauze™ (4.5 in. x 4.1 yd)
- 1x Polycarbonate Eye Shield
- 1x Combat Casualty Card
- 1x Permanent Marker



Rx ONLY

S.T.O.R.M. OPERATOR IFAK

ITEM #	COLOR	NSN
80-0133	Black	----
80-0135	Coyote	----
80-0134	OD Green	----

* Explore other IFAK's at www.NARescue.com

Learn more about this and other Products with a Mission® at www.NARescue.com

