

MARITIME GUARDIAN IFAK (MGIK)



Designed To Solve Problems, Not Create Them

The trim Maritime Guardian IFAK (MGIK) from North American Rescue® can be tucked into small spaces or carried about with ease until it is required for service. Built to the strict requirements of the United States Coast Guard, this improved first aid kit provides for safe first aid until a tactical healthcare professional arrives.

The Maritime Guardian IFAK is a reliable solution. Made of the highest quality materials, our modular system can be worn by the operator as a drop-down leg rig, hooked directly to the duty belt, or attached to an assault vest or body armor with MOLLE-style connectors.

KITS INCLUDE:



KIT CONTENTS

- 1x C-A-T® (Combat Application Tourniquet®) - Black
- 1x QuikClot® Combat Gauze®, Z-Fold Hemostatic Dressing - 3 in. x 4 yd.
- 1x NAR® S-Rolled Gauze - 4.5 in. x 4.1 yd.
- 1x ETD® Emergency Trauma Dressing - 6 in.
- 1x NAR® Nasopharyngeal Airway with Lube - 28F (117 MM)
- 1x HyFin® Vent Chest Seal (Twin Pack)
- 1x Enhanced ARS® for Needle Decompression - 10 gauge x 3.25 in.
- 1x Polycarbonate (PES) Eye Shield
- 1x GeckoGrip® Rescue Tape (6 pack)
- 2 pr Bear Claw® Nitrile Trauma Gloves - Large
- 1x Combat Casualty Card
- 1x Permanent Marker - Large



DIMENSIONS

- Case Closed: H 8.25 in. x W 71 in. x D 5.5 in.
- Weight: 2 lb 10 oz

FEATURES/SPECIFICATIONS

- Durable 500D Nylon
- Sand & water resistant, self-repairing zipper with silent 550 cord pulls
- MOLLE-style Attachment System Compatible
- Removable Belt and Leg Straps
- Leg Strap with Quick-Grip™
- Sleek Design with no Front Pocket

R_X ONLY

MARITIME GUARDIAN IFAK

| ITEM # | COLOR | NSN |
|---------|--------|------------------|
| 80-0149 | Black | 6545-01-572-7248 |
| 80-0158 | Coyote | 6545-01-587-8513 |

Learn more about this and other Products with a Mission® at www.NARescue.com



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277