

# NAR® ELITE™ - XL RESPONDER LITTER

## Evacuation Litter for Incident Tactical Extraction

North American Rescue's ELITE™ XL Responder Litter is a larger, cost-effective emergency carry litter designed for rapid deployment during casualty movement and is the ideal disposable litter for mass casualty incidents, disasters, active shooter events, as well as for everyday patient movements. This multi-use, water-resistant, disposable evacuation platform is 12 in. wider and 6 in. longer than our standard ELITE™ and features 14 hand-holds (six on each side and one larger handle at each end) plus two built-in patient retention straps for easier and safer patient handling. The ELITE™ XL is rated to 1,000 lb payload, offering more room, easier handling and overall safer patient transport.

Check out the entire series of NAR ELITE™ Litter products for the right features based on your requirements.



### FEATURES/SPECIFICATIONS

- Cost effective emergency carry litter
- Features (14) fourteen large carry handles with a load capacity tested to 1000 lb
- Includes (2) two built-in patient retention straps for easier and safer patient handling
- Latex free, water-resistant

### DIMENSIONS

- Litter Bed: L 80 in. x W 40 in.
- Packaged: L 11 in. x W 7.75 in. x D 3 in.
- Weight: 2.2 lb
- Payload: 1,000 lb



Learn more about this and other Products with a Mission® at [www.NARescue.com](http://www.NARescue.com)

### NAR® ELITE™ - XL RESPONDER LITTER

| ITEM #  | NSN  |
|---------|------|
| 60-0065 | ---- |



**NORTH AMERICAN RESCUE®**  
[www.NARescue.com](http://www.NARescue.com) • 888.689.6277