

# NAR® ELITE - LITE RESPONDER LITTER

Learn more about this and other Products with a Mission® at [www.NARescue.com](http://www.NARescue.com)



Expanding on the success and feedback from the ELITE™ and ELITE™ - XL Responder Litters, North American Rescue's ELITE™ (Evacuation Litter for Incident Tactical Extraction) - LITE Responder Litter has added capability and reduced weight and size to be able to fit into a cargo pocket. The ELITE™ - LITE is an emergency carry litter that is designed for rapid deployment for the movement of casualties for mass casualty incidents, disasters, active shooter events, as well as for everyday patient movements. Constructed of a lightweight nylon material with heavy-duty nylon webbing, the NAR ELITE™ - LITE Responder Litter is rated to 750 lb and is latex free. It features ten (10) carry handles with four (4) on each side and one (1) larger handle at either end of the litter to provide multiple carry points.

Check out the entire series of NAR ELITE™ Litter products for the right features based on your requirements.

## FEATURES/SPECIFICATIONS

- Cost effective emergency carry litter
- Features ten large carry handles with a load capacity tested to 750 lb. (340 kg)
- Made of lightweight, heavy-duty nylon
- Vacuum-packed to conserve space and fit in a cargo pocket
- Latex free, water-resistant
- Grey nylon webbing & stitching

## DIMENSIONS

- Litter Bed: L 74 in. x W 28 in.
- Packaged: L 7 in. x W 6 in. x D 17.5 in.
- Weight: 13 oz
- Load rating: 750 lb



## NAR ELITE - LITE RESPONDER LITTER

ITEM #

60-0071

\* Explore other soft litters and patient evacuation solutions at [www.NARescue.com](http://www.NARescue.com)



**NORTH AMERICAN RESCUE®**

[www.NARescue.com](http://www.NARescue.com) • 888.689.6277