

RESCUE TASK FORCE CHEST POUCH MEDICAL KIT

Our Rescue Task Force Chest Pouch is designed to add medical equipment to the MKIII RTF Tactical Responder Ballistic Vest. Made of 500D CORDURA® Solution-Dyed Nylon, this rugged chest pouch allows quick and easy access to the contents (non-tactical) from any ballistic vest with a MOLLE-style attachment system. The large, ambidextrous, premium-grade trauma shears are quickly accessible from the outside of the pouch and tethered to a retractable leash. It can be secured to the front of the body armor using a MOLLE style attachment system.

Rescue Task Force Chest Pouch Medical Kits are also compatible with NAR® Quick Detach Panels.



all KITS INCLUDE:



KIT CONTENTS

- 2 x C-A-T® (Combat Application Tourniquet®)
- 2 x NAR® Wound Packing Gauze (3 in. x 5 yd)
- 2 x Responder Emergency Trauma Dressing (ETD®) - 6 in.
- 2 x NAR® Nasopharyngeal Airway with Lube - 28F(117 mm)
- 2 x HyFin® Vent Chest Seal Twin Pack
- 3pr Blue Responder Gloves - Large
- 4 x Combat Casualty Card
- 1 x Scissor Leash™
- 1 x Trauma Shears (7.25 in.)
- 2 x Permanent Marker
- 1 x Luminous MED ID Patch
- 1 x Bag (NAR-4™ Chest Pouch)

Learn more about this and other Products with a Mission® at www.NARescue.com



FEATURES/SPECIFICATIONS

- Designed to add medical equipment to the MKIII RTF Tactical Responder Ballistic Vest
- Also compatible with other ballistic vests with MOLLE-style attachment system
- Configurable based on individual skill set and mission requirements
- Rapid access and organization is enhanced by elastic storage loops that keep items in place until needed
- Sand and water resistant zippers with silent pulls
- PALS/MOLLE-style attachment system
- Large, ambidextrous Trauma Shear configuration

CHEST POUCH MED KIT

ITEM #	COLOR	NSN
80-0406	Red	----
80-0407	Black	----
80-1109	COY	----

DIMENSIONS

- H 9.25 in. x W 8.50 Dimensions: L 7 in. x W 9 in. x D 3.5 in.
- Weight: 3.35 lb (3 lb 5.6 oz)



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277