

IFAK™ SHOULDER KIT

The Shoulder IFAK™ by North American Rescue is designed to equip the first responder with Basic or Advanced Life Support medical gear needed to address injuries common in the public safety environment. This hands-free, shoulder carry First Aid Kit is configured to address the leading causes of preventable death due to penetrating injury and blast trauma. The Shoulder IFAK™ is packaged in a rugged nylon high visibility red bag that is sand & water resistant.



FEATURES/SPECIFICATIONS

- Essential life-saving tools to address injuries caused by penetrating or blast trauma
- Compact hands-free individual first aid kit providing Basic or Advanced Life Support medical gear
- Easily identifiable rugged nylon bag with shoulder strap

IDEAL FOR:

- Individual Sports
- Backpacking
- Hiking
- Hunting
- Camping
- Vehicles
- Adventure Travel
- Surfing
- Motocross / ATV
- Mountain Biking
- Skiing
- Snowboarding

Learn more about this and other Products with a Mission® at www.NARescue.com



IFAK™ SHOULDER KIT

ITEM #	COLOR	OPTION	NSN
80-1124	Red	Basic	----
80-1125	Red	With Combat Gauze®	----
80-1126	Red	Advanced	----

DIMENSIONS

- Bag: H 9 in. x W 8.25 in. x D 2 in.
- Strap: L 17 in. x W 1 in.
- Weight: 1.45 lb




NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

IFAK™ SHOULDER KIT

all KITS INCLUDE:



+plus, by kit type:



KIT CONTENTS

- 1 x C-A-T® (Combat Application Tourniquet®) (Orange)
- 1 x NAR® Wound Packing Gauze (3 in. x 5 yd.)
- 1 x Responder Emergency Trauma Dressing (6 in. / Small)
- 1 x HyFin® Vent Compact Chest Seal
- 1 x Combat Gauze® LE with Z-Fold Hemostatic Bandage
- 2 x pair Blue Responder Gloves (Large)
- 1 x Responder Trauma Shears (Large 7.25 in.) (Black/Silver)
- 1 x Shoulder Sling IFAK Bag (Red)

BASIC KIT ALSO INCLUDES:

- 1 x NAR® Wound Packing Gauze (3 in. x 5 yd.)

BASIC WITH COMBAT GAUZE® KIT ALSO INCLUDES:

- 1 x NAR® Wound Packing Gauze (3 in. x 5 yd.)
- 1 x Combat Gauze LE with Z-Fold Hemostatic Bandage

ADVANCED KIT ALSO INCLUDES:

- 1 x NAR® Wound Packing Gauze (3 in. x 5 yd.)
- 1 x ARS® Needle Decompression Kit (10 GA / 3.25 in.)
- 1 x Combat Gauze LE with Z-Fold Hemostatic Bandage
- 1 x Pre-lubed NAR® Nasopharyngeal Airway Kit - 28F

Learn more about this and other Products with a Mission® at www.NARescue.com

