

Learn more about this and other Products with a Mission® at www.NARescue.com

K-9 Handler ROO® M-FAK

Mini First Aid Kit



One of our most compact First Aid Kits empowering K9 handlers with rapid one handed access to your tourniquet without ever having to open the full kit!

The K9 Handler ROO® M-FAK, one of North American Rescue®'s most compact Mini First Aid Kit platforms, was designed with versatile, multi-use options for delivering immediate critical care in the treatment of penetrating, blast or other traumatic injuries for K9s and their handlers in the line of duty. Despite its reduced size, this kit comes fully loaded with the critical medical equipment to address K9 emergencies from major bleeding to bloat or heat injury. Treatment pack out options include Basic and Basic with Combat Gauze® LE.

The K9 Handler ROO® M-FAK's super compact, 500 Denier nylon bag features a unique front tourniquet pouch that easily opens with the pull of a tab for rapid one handed access to your tourniquet without the need to open the full kit. Multiple elastic loops on the interior allow for secure gear organization, while it's rugged exterior includes a MOLLE style attachment for mounting in the standard vertical position as well as a 3 in. internal sleeve that accommodates horizontal attachment to your duty belt. Also included is a 1 in. PVC Med Cross Patch for convenient identification.




FEATURES/SPECIFICATIONS

- One of NAR®'s most compact, versatile Mini First Aid Kits for point of injury care for canines and/or humans
- Features rapid, one-handed access to your tourniquet from an easy-to-open front pouch allowing the C-A-T® to be used as a true one-handed tourniquet if needed
- Contains equipment for treating penetrating, blast or other traumatic injuries in the line of duty
- Designed to address K9 emergencies from major bleeding to bloat or heat injury
- Super compact, rugged nylon platform that allows attachment both vertically (MOLLE backing) and horizontally (3 in. belt loop)
- Multiple elastic loops for secure gear organization
- Includes 1 in. PVC Med Cross ID Patch
- Compatible with NAR® Quick Detach Panels (sold separately)

K-9 HANDLER ROO® M-FAK

ITEM #	COLOR	OPTION
80-1067	Black	Basic
80-1068	Coyote	Basic
80-1069	OD Green	Basic
80-1070	Ranger Green	Basic
80-1071	Black	With Combat Gauze®
80-1072	Coyote	With Combat Gauze®
80-1073	OD Green	With Combat Gauze®
80-1074	Ranger Green	With Combat Gauze®
80-1424	Black	With Celox™ Gauze
80-1425	Coyote	With Celox™ Gauze
80-1426	OD Green	With Celox™ Gauze
80-1427	Ranger Green	With Celox™ Gauze
80-1432	Black	Advanced with Combat Gauze®
80-1433	Coyote	Advanced with Combat Gauze®
80-1434	OD Green	Advanced with Combat Gauze®
80-1435	Ranger Green	Advanced with Combat Gauze®
80-1428	Black	Advanced with Celox™ Gauze
80-1429	Coyote	Advanced with Celox™ Gauze
80-1430	OD Green	Advanced with Celox™ Gauze
80-1431	Ranger Green	Advanced with Celox™ Gauze

Nylon Colors:



DIMENSIONS

- L 7 in. x W 3.75 in. x D 4 in.
- Weight: 13 oz

K-9 Handler ROO® M-FAK

Mini First Aid Kit

all KITS INCLUDE:



+plus, BASIC With:



+plus, ADVANCED With:



K-9 HANDLER ROO® M-FAK

ITEM #	COLOR	OPTION
80-1067	Black	Basic
80-1068	Coyote	Basic
80-1069	OD Green	Basic
80-1070	Ranger Green	Basic
80-1071	Black	With Combat Gauze®
80-1072	Coyote	With Combat Gauze®
80-1073	OD Green	With Combat Gauze®
80-1074	Ranger Green	With Combat Gauze®
80-1424	Black	With Celox™ Gauze
80-1425	Coyote	With Celox™ Gauze
80-1426	OD Green	With Celox™ Gauze
80-1427	Ranger Green	With Celox™ Gauze
80-1432	Black	Advanced with Combat Gauze®
80-1433	Coyote	Advanced with Combat Gauze®
80-1434	OD Green	Advanced with Combat Gauze®
80-1435	Ranger Green	Advanced with Combat Gauze®
80-1428	Black	Advanced with Celox™ Gauze
80-1429	Coyote	Advanced with Celox™ Gauze
80-1430	OD Green	Advanced with Celox™ Gauze
80-1431	Ranger Green	Advanced with Celox™ Gauze



Learn more about this and other Products with a Mission® at www.NARescue.com

