

Learn more about this and other Products with a Mission® at www.NARescue.com



M-FAK LCL

MINI FIRST AID KIT

NAR®'s most compact, versatile Mini First Aid Kit, now available in laser cut nylon!

North American Rescue's Mini First Aid Kit (M-FAK) - now featuring Laser Cut Laminate (LCL) MOLLE - was designed to be the most compact, multi-use IFAK for delivering immediate critical care for penetrating, blast or other traumatic injuries. At half the size of our popular T.O.R.K. platform, the M-FAK is currently the smallest platform available in our product offering. Despite its reduced size, this kit comes fully loaded with the critical medical equipment most requested by First Responders operating in the line of duty. Treatment pack out options include Basic, Basic with Combat Gauze® LE and Advanced with Combat Gauze® LE versions.

This Individual First Aid Kit's super compact, rugged nylon bag opens via 550 cord zipper pulls from the side in a clam shell configuration utilizing (2) two main sleeves that open on both ends for easy access and features multiple elastic loops for secure gear organization. The bag's exterior includes a MOLLE backing for mounting in the standard vertical position as well as a 3 inch internal sleeve for accommodating horizontal attachment to your duty belt.

The vertical mount can be set to open left-to-right or right-to-left based on shooter preference, while the horizontal mount on a belt allows opening directly to your C-A-T® tourniquet. Also included are a loop patch on front of the bag for custom user labeling and (3) three MOLLE strips on the front, and (2) two MOLLE strips on the spine's edge to attach additional gear such as TQ/holders, tactical lights, or knives (sold separately) if needed.

FEATURES/SPECIFICATIONS

- NAR's most compact, versatile Mini First Aid Kit - now with laser cut laminated MOLLE
- Contains the First Responder's most requested critical point-of-wounding medical equipment for treating penetrating, blast or other traumatic injuries in the line of duty
- The LCL IFAKS are double spaced for the MOLLE attachment points on the rear of each IFAK for maximum adjustability and super secure mounting to any platform
- Half the size of our popular T.O.R.K platform

DIMENSIONS

- L 6 in. x W 3 in. x D 3 in.
- Weight: 13 oz

M-FAK LCL

| ITEM # | COLOR | OPTION |
|---------|--------------|-----------------------------|
| 80-1667 | Black | Basic |
| 80-1668 | Coyote | Basic |
| 80-1669 | OCP | Basic |
| 80-1670 | OD Green | Basic |
| 80-1671 | Red | Basic |
| 80-1672 | Ranger Green | Basic |
| 80-1673 | Black | Basic with Combat Gauze® |
| 80-1674 | Coyote | Basic with Combat Gauze® |
| 80-1675 | OCP | Basic with Combat Gauze® |
| 80-1676 | OD Green | Basic with Combat Gauze® |
| 80-1677 | Red | Basic with Combat Gauze® |
| 80-1678 | Ranger Green | Basic with Combat Gauze® |
| 80-1679 | Black | Advanced with Combat Gauze® |
| 80-1680 | Coyote | Advanced with Combat Gauze® |
| 80-1681 | OCP | Advanced with Combat Gauze® |
| 80-1682 | OD Green | Advanced with Combat Gauze® |
| 80-1683 | Red | Advanced with Combat Gauze® |
| 80-1684 | Ranger Green | Advanced with Combat Gauze® |



Nylon Colors:



M-FAK LCL

MINI FIRST AID KIT

all KITS INCLUDE:



+plus, by kit type:



Learn more about this and other Products with a Mission® at www.NARescue.com



KIT CONTENTS

- 1x C-A-T® (Combat Application Tourniquet®) - Black
- 1x NAR® Wound Packing Gauze - 3 in. x 5 yd.
- 1x Hyfin® Compact (Vented) Twin Pack Chest Seal Dressing
- 1x Flat Emergency Trauma Dressing - 4 in.
- 1 pr Black Responder Gloves - Large