

BLOOD ARMOR DUO

The Blood Armor Duo, co-developed by North American Rescue® and MAXQ, is designed and validated for transporting refrigerated blood products for emergency pre-hospital transfusions. The container is capable of holding 2 units of Whole Blood or 3 units of Packed Red Blood Cells (PRBCs) - units are not included. The container is developed with packing configurations to allow the end user (military field medics, dismantled units, and civilian emergency medical services) the ability to design and support pre-hospital and on-the-go transfusion program without the traditional cooler limitations. The container is designed with a standard universal pack-out that is meant for most standard global operating conditions (Summer, Winter, Extreme Summer, and Controlled ambient). Additionally, configuration supports transport temperature stability in arctic conditions.

FEATURES/SPECIFICATIONS

- Vacuum Insulated panels are protected by the cooler's Blood Armor Duo's internal Cordura layer, adding extra durability and reducing the risk of a puncture
- Designed and validated for transporting refrigerated blood products for emergency pre-hospital transfusions
- Lightweight, designed with a standard universal pack-out meant for most standard global operation conditions
- Capability of 2 units Whole Blood or 3 units Packed Red Blood Cells (PRBCs). Units not included
- Arctic and hot weather tested, Proven protection in -33°C (-27.4°F) to 40°C (104°F)

BLOOD NOT INCLUDED



Learn more about this and other Products with a Mission® at www.NARescue.com



TEMPERATURE HOLD PERFORMANCE (1 - 10°C / 33.8 - 50°F) AT VARIOUS OPERATING CONDITIONS

Qualification Ambient	Capacity	Pack-Out	Qualification Duration
Summer Ambient (32.5°C / 90.5°F)	2 WB / 3 RBC	Universal	14 hours
Winter Ambient (ISTA 7D Winter)	2 WB / 3 RBC	Universal	14 hours
Controlled Room Temperature	2 WB / 3 RBC	Universal	24 hours
Extreme Summer Ambient (40°C / 104°F)	2 WB / 3 RBC	Universal	10.5 hours
Arctic Ambient (MIL-STD 810G Basic Cold)	2 WB / 3 RBC	Arctic	11 hours

BLOOD ARMOR DUO

ITEM #

35-0076

* Explore other blood cooler products at www.NARescue.com

DIMENSIONS

- External Dimensions: L 8 in. x W 5.5 in. x D 5.75 in.
- Internal Dimensions: L 6.75 in. x W 4.25 in. x D 4 in.
- Packaged Weight: 1.82 lb.



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277