

EAGLE IFAK

One of NAR's smallest IFAKs, complete with essential life-saving equipment for point-of-wounding and self-aid/buddy-aid care to address penetrating, blast or other traumatic injuries.

The Eagle IFAK is a compact, versatile individual first aid kit that offers both MOLLE/PALS style connectors to attach to the assault vest/body armor and an upper or lower mounted 2.25" belt loop attachment to secure the kit in the high or low position on your duty/pistol belt, all in one carrying platform. Sized as one of the smallest IFAKs in our product line, this innovative kit contains all of the essential equipment for point-of-wounding and self-aid/buddy-aid care to address penetrating, blast or other traumatic injuries. Featuring a Quick Pull Tab, the carrying bag easily opens to reveal the removable kit with just one hand. The kit has a unique, detachable shock-cord inner panel that can be adjusted to carry different components and provides easy access to your IFAK contents by simply pulling the color-coded tab to remove the panel from the carrying bag. The compact size reduces the space requirement needed for an IFAK on your vest or belt.



EAGLE IFAK

ITEM #	COLOR	OPTION	NSN
85-0177	Black	Basic	----
85-0178	Coyote	Basic	----
85-0180	OD Green	Basic	----
80-1491	RGR Green	Basic	----
80-1490	Multicam	Basic	----
85-0915	Black	Basic with Combat Gauze®	----
85-0982	Coyote	Basic with Combat Gauze®	----
85-1044	OD Green	Basic with Combat Gauze®	----
80-1492	RGR Green	Basic with Combat Gauze®	----
85-0191	Multicam	Basic with Combat Gauze®	----
85-0189*	Black	Advanced	----
85-0190*	Coyote	Advanced	----
85-0192*	OD Green	Advanced	----
80-1496*	RGR Green	Advanced	----
80-1495*	Multicam	Advanced	----

* Adv. kits are Rx ONLY

DIMENSIONS

- Packaged: H 8 in. x W 4 in. x D 3 in.
- Weight: 1 lb 3 oz

FEATURES/SPECIFICATIONS

- Can connect to the assault vest/body armor or to a duty/pistol belt
- Contains the essential equipment needed for self-aid/buddy-aid to address the top leading causes of potentially preventable death in a tactical environment
- Quick Pull Tab easily opens the kit with just one hand
- Has a unique, detachable shock-cord inner panel that can be adjusted to your internal components
- Inner panel is easily removable for easy access to your kit
- Compatible with NAR® Quick Detach Panels (sold separately)

Learn more about this and other Products with a Mission® at www.NARescue.com



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

EAGLE IFAK

all KITS INCLUDE:




KIT CONTENTS

- 1x C-A-T® - Black (Combat Application Tourniquet®)
- 1x 4 in. Flat ETD®
- 1x HyFin® Vent Chest Seal Twin Pack
- 1 pr Bear Claw® Nitrile Trauma Gloves, Lg.

+plus, by kit type:

- 1x NAR® Wound Packing Gauze - **- OR -**
- 1x Z-Fold Combat Gauze® - **- OR -**
- 1x Z-Fold Combat Gauze®
- 1x NAR® Nasopharyngeal Airway with Lube - 28F (117 mm)
- 1x ARS® Needle Decompression Kit (14 G x 3.25 in.)

+plus, by kit type:



EAGLE IFAK

ITEM #	COLOR	OPTION	NSN
85-0177	Black	Basic	----
85-0178	Coyote	Basic	----
85-0180	OD Green	Basic	----
80-1491	RGR Green	Basic	----
80-1490	Multicam	Basic	----
85-0915	Black	Basic with Combat Gauze®	----
85-0982	Coyote	Basic with Combat Gauze®	----
85-1044	OD Green	Basic with Combat Gauze®	----
80-1492	RGR Green	Basic with Combat Gauze®	----
85-0191	Multicam	Basic with Combat Gauze®	----
85-0189*	Black	Advanced	----
85-0190*	Coyote	Advanced	----
85-0192*	OD Green	Advanced	----
80-1496*	RGR Green	Advanced	----
80-1495*	Multicam	Advanced	----

* Adv. kits are Rx ONLY

Learn more about this and other Products with a Mission® at www.NARescue.com

