

RESCUE TASK FORCE VEST KIT^{with MED Kit}

REV2018-09-19

Developed for Rescue Personnel Working in the Direct and Indirect Threat Areas of Active Shooter and Mass Casualty Incidents...

The Rescue Task Force Vest Kit combines state-of-the-art ballistic personal protection along with life-saving, battlefield proven medical equipment from North American Rescue. This ideal solution was designed for EMS, Fire-Rescue and Law Enforcement personnel performing medical operations in response to Active Shooter and Mass Casualty Incidents (AS/MCI).

Developed for rescue personnel working in the Direct and Indirect Threat areas of AS/MCI's, the Rescue Task Force Vest Kit provides rescuers with a personal protective ballistic vest combined with the essential medical equipment to provide immediate point-of-wounding care to injured casualties in accordance to Tactical Emergency Casualty Care (TECC) guidelines.

The ballistic protective vest provides rescuers with one of the largest NIJ 0101.06 Level IIIA front/back panel soft armor protection systems available in the industry: Front/back panel protection area: 187 sq. in. (W 11 in. x H 17 in.). The vest comes with built-in front and back hard armor plate pouches to allow the option to upgrade up to various sizes of Level IV hard plate protection (plates sold separately). The ballistic soft armor is comprised of a unique, multi-hit capable, patent-pending assembly of ballistic materials that capture projectiles and disperse the energy over the entire surface of the panel. Front and back ballistic protection has the capability to be upgraded with the addition of Level IV hard plates (sold separately).

This "one-size-fits-most" soft armor ballistic vest easily adjusts at the shoulder and waist to fit different size rescuers. The front & back panels and shoulder straps have MOLLE capability that allows for the attachment of additional medical or extraction equipment. Designed to address the top leading causes of potentially preventable death in a tactical environment, the medical equipment can also be modified to meet your agency's needs as required by your protocols.



Front/back panel protection area: 187 sq. in. (W 11 in. x H 17 in.)

RED
BLK
YEL
COY
BLU

Rescue Task Force Vest Kit with Medical Kit

Color	Item #	NSN#
RED	80-0383	
YEL	80-0384	
BLK	80-0363	
COY	80-0364	
BLU	80-0447	



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS
WITH A
MISSION®

RESCUE TASK FORCE VEST KIT *with MED Kit*

REV2018-09-19

Features:

- NIJ 06 Level IIIA Soft Armor Ballistic Protection:
- Front and back hard armor plate pockets able to receive multiple size plates
- Mil-Spec industrial reinforced webbing
- Front and back MOLLE modular attachment system (PALS compatible)
- Durable 1,000-denier Cordura nylon external carrier construction
- Adjustable hook and loop side straps expandable up to 62 in.
- Rescue Handle – 1,200 lb tensile strength reinforced strap
- One of the largest soft armor coverage areas in the industry: Front/back panel protection area: 187 sq. in. (W 11 in. x H 17 in.)
- Adjustable hook and loop padded shoulder straps with MOLLE plus two fixed "D" rings
- 5-year manufacturer's warranty
- Blast mitigation configuration
- Comes complete with medical equipment (NAR-4 Chest Pouch with modified equipment list and additional C-A-T® tourniquets /holders)
- Made in the USA.

Dimensions:

- Vest with Medical Kit (stored): H 19 in. x W 13.5 in. x D 5 in.
- Soft Armor: H 17 in. x W 11 in.
- ID Placards: H 2 in. x W 11 in.
- Maximum Waist Size: 62 in.
- Total Weight: 9 lb

RESCUE TASK FORCE VEST MEDICAL KIT CONTENTS:



*Front & Back
"Rescue Task Force"
Placards* are included*



**Custom Agency ID Placards also available for purchase. Contact Customer Service for ordering information.*



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS

WITH A

MISSION®