

CUATTRO TM RUGGED ACQUISITION TABLET

A Smarter direction in Digital Radiography... Better. Much better.

The Cuattro Tactical Medical (TM) Rugged Acquisition Tablet is an 13.3 in. Sunlight Readable Rugged Tablet with CUATTRO Medical Software. It is designed to work with other Cuattro TM™ products. Tablet only. Other Cuattro TM products sold separately.

Cuattro provides a comprehensive 1-year warranty that includes full coverage against defects in materials and workmanship, integrated software support, and access to 24/7/365 U.S.-based Customer Support with Co-Pilot remote assistance for real-time diagnostics, evaluation, resolution of any warranty or service needs. Remote Hardware Monitoring, Advanced Loaner Program (detectors and workstations) and a 99% uptime guarantee.



FEATURES/SPECIFICATIONS

- Compatible with other Cuattro Tactical Medical (TM) products (sold separately)
- 1-Year Warranty included
- Processor: Intel 13th Gen i7
- Storage: 512 GB SSD
- Installed RAM: 16 GB
- Operating System: Windows 11
- Screen: Capacitive Touch LCD
- Additional Information: TPM 2.0 and TAA Compliant with Accessories, Power Supply. Includes Cuattro Tactical software.

DIMENSIONS

- Dimensions: L 10.8 in. x W 12.8 in. x D 0.8 in.
- Weight: 3.79 lbs.

Learn more about this and other Products with a Mission® at www.NARescue.com



CUATTRO TM RUGGED ACQUISITION TABLET

ITEM #

20-0126

* Explore other Cuattro products at www.NARescue.com