

RESCUE TASK FORCE VEST KIT *with Side Armor*

REV2018-09-19

State-of-the-Art Ballistic Personal Protection along with Life-Saving, Proven Medical Equipment...

The Rescue Task Force Vest Kit combines state-of-the-art ballistic personal protection along with life-saving, battlefield proven medical equipment from North American Rescue. This ideal solution was designed for EMS, Fire-Rescue and Law Enforcement personnel performing medical operations in response to Active Shooter and Mass Casualty Incidents (AS/MCI).

Developed for rescue personnel working in the Direct and Indirect Threat areas of AS/MCI's, the Rescue Task Force Vest Kit provides rescuers with a personal protective ballistic vest combined with the essential medical equipment to provide immediate point-of-wounding care to injured casualties in accordance to Tactical Emergency Casualty Care (TECC) guidelines.

The ballistic protective vest provides rescuers with one of the largest NIJ 0101.06 Level IIIA front/back panel soft armor protection systems available in the industry: Front/back panel protection area: 187 sq. in. (W 11 in. x H 17 in.). The vest comes with built-in front and back hard armor plate pouches to allow the option to upgrade up to various sizes of Level IV hard plate protection. The ballistic soft armor is comprised of a unique, multi-hit capable, patent-pending assembly of ballistic materials that capture projectiles and disperse the energy over the entire surface of the panel. The kit also includes the optional cummerbund-style side Level IIIA soft armor (coverage H 7 in. x L 62 in.) that has the capability of adding additional lateral protection via side plate pouches (sold separately) to accommodate hard side plates (not included). This full configuration offers Level IIIA soft armor front, back and side ballistic protection that has the capability to be upgraded with Level IV hard plates (not included).

This "one-size-fits-most" soft armor ballistic vest easily adjusts at the shoulder and waist to fit different size rescuers. The front and back panels, shoulder straps and optional cummerbund side armor have MOLLE capability that allows for the attachment of additional medical or extraction equipment. Designed to address the top leading causes of potentially preventable death in a tactical environment, the medical equipment can also be modified to meet your agency's needs as required by your protocols. Armor and medical products can also be purchased separately.



△ **Coyote**
(Coyote Medical Pouch,
TQ Holders & ID Placards)



△ **Blue**
(Black Medical Pouch
TQ Holders & ID Placards)



△ **Black**
(Black Medical Pouch, TQ Holders & ID Placards)



△ **Red**
(Red Medical Pouch,
TQ Holders & ID
Placards)

Rescue Task Force Vest Kit with Side Armor/Medical Kit

Color	Item #	NSN#
RED	80-0385	
YEL	80-0386	
BLK	80-0365	
COY	80-0366	
BLU	80-0448	



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277



NOTE: Tactical Responder Side Armor Pouches
are available for separate purchase

RESCUE TASK FORCE VEST KIT *with Side Armor*

REV2018-09-19

Yellow ▽

(Black Medical Pouch, TQ Holders & ID Placards)



1
Rescue Responder
Soft Side Armor*

2
C-A-T® Holders
(Left & Right)

*NOTE: Tactical Responder Side Armor
Pouches are available for separate purchase

RESCUE TASK FORCE VEST MEDICAL KIT CONTENTS:



Features:

- NIJ 06 Level IIIA Soft Armor Ballistic Protection
- Front and back hard armor plate pockets able to receive multiple size plates
- Includes cummerbund-style Side Armor Protection (coverage H 7 in. x L 62 in.) with the capability of additional plate pouches (sold separately)
- Mil-Spec industrial reinforced webbing
- Front and back MOLLE modular attachment system (PALS compatible)
- Durable 1,000-denier Cordura nylon external carrier construction
- Adjustable hook and loop side straps – expandable up to 62 in.
- Rescue Handle – 1,200 lb tensile strength reinforced strap
- One of the largest soft armor coverage areas in the industry – Front/back panel protection area: 187 sq. in. (W 11 in. x H 17 in.)

- Adjustable hook and loop padded shoulder straps with MOLLE plus two fixed “D” rings
- 5-year manufacturer’s warranty
- Blast mitigation configuration
- Comes complete with medical equipment (NAR-4 Chest Pouch with modified equipment list and additional C-A-T® tourniquets /holders)
- Made in the USA

Dimensions:

- Vest with Medical Kit & Side Armor (stored): H 19 in. x W 13.5 in. x D 5 in.
- Front & Back Soft Armor panels: H 17 in. x W 11 in.
- Cummerbund coverage area: H 7 in. x L 62 in.
- ID Placards: H 2 in. x W 11 in.
- Maximum Waist Size: 62 in.
- Total Weight: 13 lb

****Custom Agency ID Placards also available for purchase.
Contact Customer Service for ordering information.**



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS
WITH A
MISSION®