

NAR-6 KIT

Low-Profile Tactical Medical Kit for Combat Medics

The NAR-6 Aid Kit is a mission-ready, low-profile **tactical medical kit** engineered to ensure military personnel and first responders maintain immediate access to critical **life-saving equipment** in austere and high-threat environments. Designed for versatility, this compact **combat medic kit** supports a wide range of mission profiles, from point-of-injury care to prolonged field operations. Constructed from rugged 500D CORDURA®, this durable **medic bag** is built to perform in demanding conditions while maintaining a streamlined, low-profile footprint. Sand- and water-resistant, self-repairing zippers with silent 550 cord pulls enable quiet operation during **tactical movements**. The interior features six modular hook-and-loop pouches in varying sizes, allowing for full **mission-specific loadout customization**. Four smaller pouches are optimized for **hemorrhage control** supplies such as tourniquets, gauze, and pressure bandages, while two larger pouches accommodate **airway management** and **IV/IO medical equipment**.

Detachable internal MOLLE panels enhance organization and provide rapid access to **essential trauma supplies**, enabling efficient "grab-and-go" capability in high-stress scenarios. Integrated mesh compartments offer additional storage for **medical accessories** and supplemental gear. An externally accessible bottom pouch delivers **immediate IFAK (Individual First Aid Kit) capability**, allowing rapid deployment of **critical bleeding control tools** without opening the main compartment.

Each removable pouch is designed for quick identification and direct use at the point of injury, supporting **efficient casualty care** and **TCCC (Tactical Combat Casualty Care) protocols**. Hook-and-loop attachment ensures secure retention while enabling rapid detachment for focused treatment of specific **life-threatening injuries**.

The slim, low-profile design minimizes bulk while maximizing operational efficiency. The kit can be secured with anti-pilferage seals to protect sensitive **medical supplies** during transport or staging. Stowable shoulder straps allow the bag to convert quickly into a backpack-style **field medical pack**, supporting multiple carry configurations based on mission requirements.

Optimized for **military medics, tactical operators, and first responders**, the NAR-6 Aid Kit delivers a scalable, durable, and highly organized solution for **combat casualty response, trauma care, and emergency medical preparedness**.

FEATURES/SPECIFICATIONS

- Durable 500D CORDURA® construction for combat-ready performance
- Low-profile, ergonomic design for streamlined tactical loadouts
- Six removable, mission-configurable internal pouches
- Optimized for IFAK setups, trauma kits, and combat medic applications
- MOLLE-compatible internal panels for modular organization
- External quick-access IFAK pouch for rapid deployment
- Silent zipper pulls for reduced noise signature
- Stowable shoulder straps for backpack-style carry
- Compatible with anti-pilferage seals for gear security



NAR-6 KIT

R_X ONLY

ITEM #	COLOR	NSN
80-1972	Black	----
80-1973	Coyote	----

* Explore the NAR-6 Kit at www.NARescue.com

DIMENSIONS


- Dimensions: L 20 in. x W 12.5 in. x D 6 in.
- Weight: 10.5 lbs.



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

NAR-6 KIT

all KITS INCLUDE:



Learn more about this and other Products with a Mission® at www.NARescue.com



KIT CONTENTS

- 5 x C-A-T® Tourniquet - Black
- 2 x C-A-T® Holder (LCL) - Black **or** Coyote
- 2 x QuikClot® Z-Fold Combat Gauze® - 3 in. x 4 yd.
- 4 x NAR® Wound Packing Gauze - 3 in. x 5 yd.
- 2 x NAR® Flat Emergency Trauma Dressing (ETD®) - 4 in.
- 1 x NAR® Emergency Trauma Dressing (ETD®) - 6 in.
- 1 x Abdominal/Stump Emergency Trauma Dressing (ETD®) - 12 in. x 12 in.
- 2 x NAR® Nasopharyngeal Airway with Lube - 28F (117mm)
- 1 x Tactical Suction Device
- 1 x Cyclone® Pocket BVM
- 1 x BAC-Pack Bougie-Aided Cricothyroidotomy Pack
- 2 x Enhanced ARS Needle - 10G
- 4 x HyFin® Vent Chest Seal Dressing - Twin Pack
- 2 x NAR® Survival Blanket
- 2 x Elastic Wrap Bandage - 4 in.
- 2 x BurnTec® Dressing - 4 in. x 4 in.
- 1 x SAM® Splint II
- 2 x Triangular Bandage - OD Green
- 1 x Gecko Grip® Flat Tape (6pk) - 2 in.
- 1 x Minor Care Bandages Kit
- 1 x Supplemental Aid Module Kit
- 2 x Kerlix™ Gauze - 4.5 in. x 4.1 yd.
- 2 x PES Eye Shield
- 2 x NAR® Saline Lock Kit
- 2 x Eye Wash Solution - 4 oz.
- 4pr Black Talon Nitrile Gloves - Large
- 1pr NAR® Trauma Shears - Small
- 1 x BLEED Nametape
- 1 x AIRWAY Nametape
- 1 x BURN Nametape
- 1 x MISC. MEDS Nametape
- 1 x Eye Nametape
- 1 x Permanent Marker - Large



NAR-6 BAG

ITEM #	COLOR	NSN
80-1792T	Black	----
80-1794T	Coyote	----

* Explore the NAR-6 Kit at www.NARescue.com