

# NAR® ELITE LITTER™ - RESPONDER RED

## Evacuation Litter for Incident Tactical Extraction

North American Rescue's ELITE Litter™ - Responder Red is a cost effective emergency carry litter designed for rapid deployment for the movement of casualties and is the ideal disposable litter for mass casualty incidents, disasters, active shooter events, as well as for everyday patient movements. Constructed of a latex free water-resistant, woven polyethylene with signature red nylon webbing & stitching, the NAR® ELITE Litter™ - Responder Red is rated for up to 500 lb. and features ten (10) carry handles with four (4) on each side and one (1) larger handle at either end of the litter to provide multiple carry points. Check out the entire series of NAR ELITE Litter™ products for the right features based on your requirements.

Learn more about this and other Products with a Mission® at [www.NARescue.com](http://www.NARescue.com)



### FEATURES/SPECIFICATIONS

- Cost effective emergency carry litter
- Features ten large carry handles with a load capacity tested to 500 lb
- Latex free, water-resistant
- Signature red nylon webbing & stitching

### DIMENSIONS

- Litter Bed: L 74 in. x W 28 in.
- Packaged: L 11 in. x W 7.75 in. x D 2 in.
- Weight: 13.3 oz



### ELITE LITTER™ - RESPONDER RED

ITEM #

60-0080

\* Explore other soft litters and patient evacuation solutions at [www.NARescue.com](http://www.NARescue.com)