

NAR® ELITE LITTER™ - RESPONDER

Evacuation Litter for Incident Tactical Extraction

North American Rescue's Evacuation Litter for Incident Tactical Extraction (ELITE) Litter - Responder is a cost effective emergency carry litter designed for rapid deployment for the movement of casualties and is the ideal disposable litter for mass casualty incidents, disasters, active shooter events, as well as for everyday patient movements. Constructed of a water-resistant, woven polyethylene that is latex free and rated for up to 500 lb., the NAR® ELITE Litter™ - Responder features ten (10) carry handles with four (4) on each side and one (1) larger handle at either end of the litter to provide multiple carry points. Check out the entire series of NAR® ELITE Litter™ products for the right features based on your requirements.



FEATURES/SPECIFICATIONS

- Cost effective emergency carry litter
- Features ten large carry handles with a load capacity tested to 500 lb
- Latex free, water-resistant

Learn more about this and other **Products with a Mission®** at www.NARescue.com



NAR® ELITE LITTER™ - RESPONDER

ITEM #	NSN
60-0061	----

DIMENSIONS

- Litter Bed: L 76 in. x W 30 in.
- Packaged: L 11 in. x W 7.75 in. x D 2 in.
- Weight: 17 oz
- Payload: 500 lb