

°MEQU °M WARMER SYSTEM - POWER PACK PLUS

A Portable Blood and IV Fluid Warming System

The °MEQU Power Pack Plus (or Power Pack+) is a rechargeable battery pack that is compatible with the °MEQU °M Warmer (sold separately). A charged Power Pack Plus can warm 2 liters of cold blood or fluids (39°C +/- 3°C at flowrates from 2-150ml/min) on a single charge of the battery.

FEATURES/SPECIFICATIONS

- Rechargeable battery pack with a recharge time of 3 hours
- Compatible with the °MEQU °M Warmer (sold separately)

OPERATING PARAMETERS

- Output temperature range: 39°C (+/- 3°C at flowrates from 2-150ml/min)
- Max fluid warming volume on one fully charged Power Pack Plus (output temperature of min. 36°C:

| | |
|-------------------------|-----------------|
| Input temperature 5°C: | Power Pack+: 2L |
| Input temperature 21°C: | Power Pack+: 4L |

- Priming volume: 3.5 ml (including tubing and luer locks)
- Maximum inlet pressure: 300 mmHg

RECHARGEABLE POWER PACK PLUS SPECIFICATIONS

- Standard recharge time (Power PackPlus): 3 hours
- Standard charger input voltage AC Charger: 100-240 VAC
- Power Pack Plus: use life 2 years or 200 charging cycles
The total capacity with the number of charge cycles:

| | |
|---------------------|-------|
| -100 charge cycles: | 91.7% |
| -200 charge cycles | 83% |

- Electrical compliance: ISO 60601-1 | EN 13718-1

°MEQU products are restricted to DoW and Federal Purchases only. Please contact Customer Service for details.

°Caution: Federal (USA) law restricts this device to sale by or on the order of a practitioner licensed to use or order use of the device or Federal (USA) law prohibits dispensing without a prescription.



°MEQU °M Warmer System sold separately

°MEQU

DIMENSIONS

- Size: W 3.54 in. x L 7.09 x H 1.38 in. (9 cm. x 18 cm. x 3.5 cm.)
- Cable length: 47.24 in. (120 cm.)
- Weight: 1.57 lbs. (710 g.)

°MEQU °M WARMER SYSTEM - POWER PACK PLUS

ITEM #

35-0049

°Explore other °MEQU products at www.NARescue.com



Learn more about this and other Products with a Mission® at www.NARescue.com



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277