

NAR HEAD IMMOBILIZER

The NAR Head Immobilizer is amazingly easy to use, compact and extremely effective. Constructed of soft, closed-cell foam with hook and loop padded straps, the NAR Head Immobilizer provides quick and effective cervical restriction. Ears and head are accessible without removing the device for ongoing patient assessment.

Learn more about this and other Products with a Mission® at www.NARescue.com



Spine board sold separately



FEATURES/SPECIFICATIONS

- Provides quick and effective cervical restriction
- Soft, closed-cell foam with hook and loop padded straps
- Ears and head are accessible without removing the device
- SEL#: 09ME-04-SPIN

NAR HEAD IMMOBILIZER

ITEM #	NSN
50-0020	6530-01-619-1777

DIMENSIONS

- H 10 in. x W 16 in. x D 5.75 in.
- Weight: 2 lb 2.8 oz