

PATIENT MONITORING SUITE

WITHOUT NAR-S TABLET

The Patient Monitoring Suite by North American Rescue® Simulation is a modular and fully integrated training system designed to simulate real-time vital sign monitoring. This system enhances realism in medical training by allowing responders to interact with accurate, lifelike monitors without requiring instructor intervention. It is compatible with human role-players or any patient simulator across all phases of care.

The North American Rescue Simulation (NAR-S) Capnograph Simulator is designed to be used as a mainstream carbon dioxide monitor that can be connected in-line with a mask or tube. Simulator displays include End-tidal CO₂, respiratory rate and CO₂ waveform.

NAR Simulation's Pulse Oximeter Simulator utilizes the widely used finger-style probe for monitoring blood oxygen saturation. Simulator display includes Heart Rate, SpO₂ and SpO₂ waveforms.

The Non-Invasive Blood Pressure (NIBP) Monitor Simulator is designed to replicate a portable NIBP monitor with onboard controls. Displays include Systolic and Diastolic Pressure (mmHg), and Heart Rate.

Controlled through NAR-S tablet/Proprietary Software via Radio Signals. **NAR-S Tablet is NOT included.**

PMS INCLUDES:



FEATURES/SPECIFICATIONS

- Allows medical responders to gather vital signs without instructor intervention
- Compatible with Role Players and Manikins
- Remotely change patient vitals via NAR-S Tablet
- Includes Simulated Non-Invasive Blood Pressure Cuff, Pulse Ox, and Capnograph
- Does NOT include NAR-S control tablet

Learn more about this and other Products with a Mission® at www.NARescue.com



SUITE CONTENTS

- 1 x Pulse Oximeter Simulator
- 1 x Capnograph Simulator
- 1 x NIBP Cuff Simulator
- Hard plastic storage and transport case
- Charging system
- 3-Year Warranty

DIMENSIONS

- Dimensions: L 19 in. x W 15 in. x D 7 in.
- Weight: 9 lbs.

PMS

ITEM #

93-0263

* Explore more simulation solutions at www.NARescue.com



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277