

RIGID GEN 7 CAT® TQ CASE with Shirt Shield

REV070218

Always have the Gen7 C-A-T® within easy reach, yet keep it and your uniform shirt protected...


The Gen7 C-A-T® RIGID TQ Case® with Shirt Shield was designed to provide personnel with the capability of carrying the life-saving C-A-T® Tourniquet on their person by easily mounting it to their duty/patrol/pistol belt, where it is needed most. The included Shirt Shield also helps to prevent the undue wearing of the uniform shirt caused by the windlass clip when the two rub together. This convenient case can be mounted either vertically or horizontally and can also be attached to a Safariland Drop Leg Shroud using the existing holes (*the case mounts with 2 screws in the middle of the case*). The Injection Molded Nylon Case provides strength and durability as it holds the Gen7 C-A-T® Tourniquet tightly in place with no straps or covers. In addition, the Shirt Shield easily attaches with two screws and can be replaced if ever needed. The tourniquet can be quickly and easily removed by grabbing and pulling the top part of the C-A-T® Tourniquet. The case also includes a hook and loop cross to identify it as medical equipment. *Accommodates the Gen7 C-A-T® Tourniquet (not included)*

Features:

- Allows you to carry a Gen7 C-A-T® Tourniquet on your duty/patrol/pistol belt so you have life-saving equipment on your person at all times
- Includes a Shirt Shield complete with mounting hardware to protect your uniform shirt from undue wear caused by the friction from the windlass clip
- Blade-Tech TekLok Belt Attachment fits up to a 2.25 in. belt
- Can be mounted either vertically or horizontally or attached to a Safariland Drop Leg Shroud
- Injection Molded Nylon Case to stand up to rugged LE/MIL usage
- Cross medical identifier patch included
- NOTE: Compatible with the C-A-T® Gen7 (sold separately)

Dimensions:

- H 5.88 in. x W 2.5 in. x D 2 in.
- Weight: 5.36 oz



Rigid Gen7 C-A-T® TQ Case with Shirt Shield

Color	Description	Item #	NSN#
BLK	Black	80-0944	
COY	Coyote	80-0945	
RGR	Ranger Green	80-0946	
BWB	Black Basket Weave	80-0962	
PBK	Polished Black	80-0963	



C-A-T® Gen7
(sold separately)



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

PRODUCTS
WITH  **MISSION®**