

# WALK®

## WARRIOR AID AND LITTER KIT

Learn more about this and other Products with a Mission® at [www.NARescue.com](http://www.NARescue.com)



Includes Talon®II 90C quad-fold litter

### The WALK® is the vehicle life-saving kit for the Department of Defense.

North American Rescue® pays attention. We listen. We learn. We act. We develop products like the Warrior Aid and Litter Kit (WALK®), our innovative and reliable response to requests from Warfighters and Tactical Healthcare Professionals to re-think the purpose of vehicle life-saving kits and re-engineer them to decrease preventable death.

Success is confirmed. The US Army selected the WALK® as its official vehicle life-saving kit. How did we do it?

First, we gathered as much information as possible to correctly define the requirement. From numerous observations and after-action reviews compiled from operations Enduring and Iraqi Freedom, we learned that a new focus was needed for vehicle first-aid kits. From Warfighters and Tactical Healthcare Professionals, we learned that the focus ought to be decreasing preventable combat death at the point of wounding.

With the task defined, we developed the product. The Warrior Aid & Litter Kit (WALK®) allows first-aid treatment and immediate evacuation of one patient in a non-traditional vehicle, plus it allows multiple other trauma patients to be treated until additional transport arrives. That's successful engineering.

Our engineers took the state-of-the-art Talon 90C (a quad-fold rigid litter with retractable handle system) and designed a package that holds not only the litter, but also a streamlined pocket-system containing self- and buddy-aid medical supplies. The WALK® gives true definition to the term aid and litter kit. The WALK® also contains our Hypothermia Prevention & Management Kit (HPMK) to help stabilize patients during transport.

Light and compact, the WALK® reduces cubic weight and, thereby, evacuation time. The WALK® interfaces with any standard or non-standard casualty evacuation asset and is ideal as part of the pre-positioning of medical supplies or cache for disaster response. The WALK® has adjustable shoulder straps for occasions when the kit is used to provide treatment and evacuation on target. Bags may also be purchased separately.

North American Rescue® is customer-centric. We don't develop solutions to impose on demands. We develop solutions to meet demands.

### FEATURES/SPECIFICATIONS

- Official US Army vehicle life-saving kit
- For treatment of multiple trauma patients and immediate evacuation of one patient
- Hypothermia Prevention & Management Kit™ helps stabilize patients during transport
- Includes Talon®II 90C quad-fold litter
- Interfaces with any standard or non-standard TACEVAC asset
- Accepts both: Talon®II 81C (Retired) & Talon®II 90C Collapsible Handle



### WALK®

ITEM #	COLOR	NSN
80-0097	Black	----
80-0089	Coyote	6545-01-587-1199
80-0059	DUC	6545-01-532-4962
80-0060	OD Green	----
80-1564	Ranger	----
80-1563	Multicam	----

### DIMENSIONS

- Case Closed: H 22.75 in x W 16.5 in x D 11.5 in
- Weight: 29 lb 15 oz

\* Bag only option is available at [www.NARescue.com](http://www.NARescue.com)




**NORTH AMERICAN RESCUE®**  
[www.NARescue.com](http://www.NARescue.com) • 888.689.6277

# WALK®

## WARRIOR AID AND LITTER KIT

### all KITS INCLUDE:



### KIT CONTENTS

- 2 x C-A-T® - Black (Combat Application Tourniquet®)
- 6 x Trauma Dressing, 6 in
- 4 x S-Rolled Gauze (4.5 in x 4.1 yd)
- 1 x ETD® Abdominal (Emergency Trauma Dressing)
- 5 pr Black Talon® Nitrile Trauma Gloves
- 2 x NAR® Nasopharyngeal Airway with Lube - 28F (117 mm)
- 1 x Casualty Equipment Bag
- 2 x HyFin® Chest Seal
- 2 x ARS® Needle Decompression Kit (14 G x 3.25 in)
- 2 x SAM® Splint II
- 1 x Trauma Shears (7.25 in.)
- 1 x Surgical Tape (2 in.)
- 1 x GeckoGrip® Multipurpose Tape
- 6 x Polycarbonate Eye Shield
- 1 x Combat Casualty Reference Card
- 2 x Combat Casualty Card (Triage)
- 1 x Aviation Panel (Recognition, Orange, 35 in x 29 in)
- 1 x Talon 90C® Collapsible Handle Litter
- 1 x Hypothermia Prevention & Management Kit™ (HPMK™)
- 4 x Tie Down Strap (Universal Litter)
- 1 x Bag (CCRK®-WALK®)



Learn more about this and other Products with a Mission® at [www.NARescue.com](http://www.NARescue.com)

