

WALK® (BAG ONLY)

WARRIOR AID AND LITTER KIT

Learn more about this and other Products with a Mission® at www.NARescue.com



The WALK® is the vehicle life-saving kit for the Department of Defense.

The WALK® Bag has stream-lined pocket system built around the capacity to carry the quad-fold 90C Talon® Litter.

Complete internal pocket system and backpack straps provide easy transport from the vehicle to the point of need. The top can be un-clipped and un-rolled for quick access to the top of the compartment.

Meets NAR's 11 point Berry Amendment nylon standard in 1000 Denier Nylon standard.



Designed to fit the Talon® II 90C quad-fold litter

FEATURES/SPECIFICATIONS

- The WALK® Bag has stream-lined pocket system built around the capacity to carry the quad-fold Talon® II 90C Litter
- The top can be un-clipped and un-rolled for quick access to the top of the compartment
- The primary compartment fully unzips down each side for complete, full-length access

WALK®

ITEM #	COLOR	NSN
60-0015	Black	NATO*
60-0013	Coyote	6530-01-504-9056
60-0020	DUC	----
60-0012	OD Green	6530-01-531-6902
60-0016	Ranger	----
60-0017	Multicam	----

* NATO NSN AVAILABLE

DIMENSIONS

- H 22.75 in. x W 16.5 in. x D 11.5 in.
- Weight (bag only): 3.6 lb



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277