

LARYNGOSCOPE SET

SPECIAL OPERATIONS

Designed To Operate Effectively Under The Most Extreme Environmental Conditions

The lightweight, compact, and highly innovative Special Operations Laryngoscope Set from North American Rescue® is designed to operate effectively under the most extreme environmental conditions. Professionals rely on it. To ensure the laryngoscopes perform as long and as well as necessary, the sets are packaged with spare batteries.

The fiber optic handle is water-proof, sand-proof, and dust-proof, and features a textured gripping surface. For convenience, the set includes single-use fiber optic blades made from an enhanced polymer for exceptional strength and rigidity. For added convenience, the handle's removable light source is designed to work with any standard (green-spec) fiber optic blade.

The complete set rolls up for easy storage and easy transport. With products from North American Rescue®, being prepared takes minimal effort.



FEATURES/SPECIFICATIONS

- Lightweight & compact
- Resistant to dust, sand and water
- Operates in severe environmental conditions
- Textured, water-proof & dust-proof handle
- 80 Lumen LED examination lamp
- Removable LED light source works with any standard fiber optic blade
- Single-use fiber optic blades of enhanced polymer for exceptional strength and rigidity

SPECIAL OPERATIONS LARYNGOSCOPE SET

ITEM #	NSN
10-0001	6515-01-450-9790

DIMENSIONS

- Open: L 8 in. x W 5.5 in. x D 2 in.
- Closed: L 8 in. x W 2 in. x D 2 in.
- Weight: 10 oz.


Learn more about this and other Products with a Mission® at www.NARescue.com



LARYNGOSCOPE SET

SPECIAL OPERATIONS

KIT INCLUDES:



KIT CONTENTS

- 1 x Laryngoscope Pocket Light™ with two AAA Batteries (Item# 10-0074)
- 1 x Laryngoscope Handle (Item# 10-0005 | NSN# 6515-01-460-4736)
- 1 x Miller #3 Laryngoscope Blade (Item# 10-0011 | NSN# 6515-01-460-4681)
- 1 x Mac #3 Laryngoscope Blade (Item# 10-0013 | NSN# 6515-01-460-4685)
- 1 x Laryngoscope Pouch (Item# 10-0006 | NSN# 6515-01-467-7624)
- 2 x spare AAA batteries
- 1 x Laryngoscope Set Booklet IFU (rev. 02/15/24)

*NOTE: To order replacement parts, please call Customer Service at 1-888-689-6277

Learn more about this and other Products with a Mission® at www.NARescue.com



SPECIAL OPERATIONS LARYNGOSCOPE SET

ITEM #	NSN
10-0001	6515-01-450-9790