

# DELTA ICE TEMPWATCH

Delta Development's Tempwatch is designed to interface with the Delta ICE 2L Temperature Control Insert (sold separately) during conditioning to monitor the embedded temperature sensors and alert the user with a unique alarm melody once the PCM pack has reached 2°C. The safe temperature range for blood storage is between 2-6°C.



**Delta ICE Temperature Control Insert (2L) Sold Separately**

Learn more about this and other Products with a Mission® at [www.NARescue.com](http://www.NARescue.com)



## FEATURES/SPECIFICATIONS

- Set and Forget - Automatically Conditions the PCM Pack, taking out any Guesswork
- Displays the Current Status of PCM Pack (Too Cold, In-Range, Too Warm, Error)
- Automatically Alerts the User once the Pack is Conditioned to 2°C
- Maximizes the Holdover Time for Each Conditioned PCM Pack
- Increases Efficiency
- Water Resistant
- Designed for Bench-Top Operation
- Interfaces with Delta ICE 2L Temperature Control Insert (sold separately)
- Operating Temperatures: -20 - 50°C  
Storage Temperatures: -20 - 70°C
- USBC Cable and AC/DC Adapter Included

## DELTA ICE TEMPWATCH

ITEM #	NSN
35-0045	----

\* Explore more blood cooling solutions at [www.NARescue.com](http://www.NARescue.com).

## DIMENSIONS

- Dimensions: L 10 in. x W 10 in. x H 2 in.
- Unit Weight 1 lb.