

DELTA ICE TEMPERATURE CONTROL INSERT

2 LITER

Optimal temperature storage in virtually any environment.

The Delta ICE's Temperature Controlled Insert is an aluminum insert featuring integrated, built-in NIST traceable sensors that promote temperature uniformity and an even thaw, eliminating hot spots. This SMART technology maintains two (2) 500ml of whole blood for 48 hours at 2-6°C and eliminates the need for glycol.



* Blood not included



Learn more about this and other Products with a Mission® at www.NARescue.com



FEATURES/SPECIFICATIONS

- For use with the Delta ICE Smart Blood Cooler (sold separately)
- Aluminum insert featuring integrated, built-in NIST traceable sensors
- Maintains two 500 ml of whole blood for 48 hours at 35.6 - 42.8°F, and eliminates the need for glycol

DIMENSIONS

- Dimensions: H 5.75 in. x W 6.25 in. x D 8.24 in.
- Weight: 5 lbs.

DELTA ICE TEMPERATURE CONTROL INSERT

ITEM #	NSN
35-0043	----