

TECHLINE TRAUMA WEARABLE NECK WOUND SIMULATOR

The Techline Trauma® Wearable Neck Wound Simulator from North American Rescue® Simulation and Techline Technologies, Inc., is a silicone prosthetic device enabling live action role players to simulate a real-world neck injury. The soft, flexible, light toned skin wraps comfortably around the neck and attaches with the sewn hook and loop connectors. This realistic injury features exposed tissue and muscle and includes an integrated bleeding system tube that allows role players to demonstrate active bleeding during scenarios. Enhancing realism creates a valuable training event that better prepares first responder trainees for real-world situations. Easily cleaned and reusable, this tool enables rapid preparation for training exercises and minimizes moulage requirements – while maximizing realism.

The Techline Trauma® Wearable Neck Wound is manufactured by Techline Trauma® and is compatible with all TOMManikin® and TAMILkin® simulators.



Learn more about this and other Products with a Mission® at www.NARescue.com



FEATURES/SPECIFICATIONS

- Hook and loop closure
- Exposed tissue and muscle
- Integrated bleeding system

DIMENSIONS

- L 6 in. x W 6 in. x H 6 in.
- Weight: 1 lbs.

WEARABLE NECK WOUND

ITEM #	NSN
93-0215	----

* Explore more Techline Trauma products at www.NARescue.com



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277