


TECHLINE TRAUMA WEARABLE LEG WOUND SIMULATOR

GSW

The Techline Trauma® Wearable Gunshot Wound Simulator features a through and through gunshot wound exposing fat and muscle tissues. Simulates bleeding from the entry and exit points. The hook and loop panel keeps the simulator secure during live training. This simulator is supplied with the bleeding system. Wrap around device with hook and loop fastener.

Manufactured by Techline Trauma®

Learn more about this and other Products with a Mission® at www.NARescue.com



FEATURES/SPECIFICATIONS

- Simulates a through and through gunshot
- Bleeds from entry and exit wounds
- Washable and reusable

WEARABLE LEG WOUND - GSW

| ITEM # | NSN |
|---------|------|
| 93-0212 | ---- |

* Explore more Techline Trauma products at www.NARescue.com

DIMENSIONS

- Packaged - L 6 in. x W 6 in. x H 6 in.
- Weight: 3 lbs



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277