

BLACK TALON® GLOVES

Ultimate Nitrile Trauma Gloves

Black Talon Ultimate Nitrile gloves are made of 100% nitrile material which is free of the allergens found in standard latex gloves. They are engineered to give our users high performance by offering both maximum protection and unsurpassed tactile sensitivity. This careful attention to detail includes textured fingertips to help with fine motor skills, beaded cuffs for quick and easy donning and a thickened palm for durability and reliability. These ambidextrous gloves are packaged 50 pairs per box (Extra large, 45 pairs per box); 10 boxes per case. Designed for single-use only.

The 25 Pack contains the same great Black Talon Ultimate Nitrile gloves individually pre-wrapped in pairs and packaged for easy, rapid deployment. These ambidextrous gloves come packaged 25 pairs per bag.



Learn more about this and other Products with a Mission® at www.NARescue.com



FEATURES

- 100% nitrile
- 25 Pack contains individually wrapped pairs for convenient rapid deployment
- Box contains 100 gloves (XL, 90 gloves)
- Textured for improved wet and dry gripping
- Longer and thicker than standard nitrile gloves
- Form fitting
- Extended, beaded cuff for quick and easy donning
- Tactical Black for reduced exposure when stealth is required
- Available in sizes S, M, L and XL
- Non-sterile

GLOVE SPECIFICATIONS:

- Length: 11.8 in.
- Cuff Thickness (mils): 4.7
- Palm Thickness (mils): 6.0
- Finger Thickness (mils): 8.0

BLACK TALON® GLOVES

ITEM #	OPTION	NSN
ZZ-0218	25 Pr (S)	----
ZZ-0217	25 Pr (M)	6515-01-672-7977
ZZ-0216	25 Pr (L)	6515-01-632-0059
ZZ-0215	25 Pr (XL)	6515-01-672-7985
70-0001	Box of 100 (S)	6515-01-515-0195
70-0002	Box of 100 (M)	6515-01-515-0200
70-0003	Box of 100 (L)	6515-01-515-0197
70-0004	Box of 90 (XL)	6515-01-515-0190

* Explore more PPE at www.NARescue.com

DIMENSIONS

- Box of 50 Pair:
- Packaged: L 11 in. x W 5 in. x D 3.75 in.
 - Weight: 2 lb 4 oz
- Pack of 25 Pair:
- Packaged: L 4.75 in. x W 6 in. x D 4 in.
 - Individual pair: L 2.5 in. x W 1 in. x D .8 in.
 - Weight: 16.1 oz

