

Learn more about this and other **Products with a Mission®** at www.NARescue.com



CLASS A & B KIT Hard Cases - ANSI

A portable solution to store & carry First Aid equipment in the workplace or at a work site

This rugged hard case is easily transported to the patient's side and protects contents from the elements, dust and insects. It closes securely using latches, utilizes a gasket to securely seal it, and includes a pressure release valve to ensure that it will easily open if moved between high and low pressure environments. The case offers two locking points that can be used to secure it, or employ tamper seals to ensure contents can be monitored for use. Inner nylon pouches with zipper closures are held in position using hook and loop. Each pouch has a hook and loop patch on the lid for easy identification of the contents. The Class B case has a lid that offers clear plastic pouches for easily viewable storage. These cases offer users the ability to customize the position and contents to suit their specific need.



Class A Case

**Case and
inner bags only**



Class B Case



CLASS A & B HARD CASES - ANSI

ITEM #	SIZE	COLOR	NSN
80-1113	Class A	RED	----
80-1871	Class A	ORG	----
80-1872	Class A	YEL	----
80-1114	Class B	RED	----
80-1873	Class B	ORG	----
80-1874	Class B	YEL	----

DIMENSIONS

Class A Case

- H 6 in. x W 13.5 in. x D 11.5 in.
- Weight: 5.89 lbs.

Class B Case

- H 6 in. x W 13.5 in. x D 11.5 in.
- Weight: 8.87 lbs.

FEATURES/SPECIFICATIONS

- Modular Type IV ANSI Container featuring a rugged storage organizer in the lid and removable nylon pouches with easy-to-read treatment labels
- Rugged hard case that protects your medical equipment from the elements
- Designed to be easily carried to the patient's side
- Features pressure equalization valve for air transport
- Can be locked or tamper sealed



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277