

Learn more about this and
other Products with a Mission®
at www.NARescue.com



NAR Quick Detach Panels

2 size options

Allows the operator the fastest access to their kit during the critical moments after wounding...

The new line of NAR Quick Detach Panels from North American Rescue allows the operator the fastest access to their kit during the critical moments after wounding. Utilizing Mil-spec hook & loop and heavy-duty metal snaps to ensure that the detachable front panel remains secure in the most austere of conditions, the Quick Detach Panel can be used with any Mil-spec standard Molle/PALS system. Three versions of the panel designed to compliment all of NAR IFAK's are available including a NEW horizontal mount option



FEATURES/SPECIFICATIONS

- Allows any Molle/PALS kit to be detachable
- Gives user faster access to Molle/PALS equipped kits
- Mil-Spec hook & loop and heavy duty snaps to ensure kit security
- Designed to be compatible with all of NAR IFAK's, including a NEW horizontal mount option in the medium version (W 4.5 x H 6.5 in.)

NAR QUICK DETACH PANEL (3 X 6 IN.)



DIMENSIONS

- W 3 in. x H 6 in. x D 0.5 in.
- Weight: 2.8 oz

ITEM #	COLOR
80-0893	Black
80-0894	Coyote

FEATURES/SPECIFICATIONS

- NAR's smallest detach panel measuring (W 3 x H 6 in.)
- Utilizes only 3 rows by 2 columns on user vest/bag, vertical mount configuration
- Designed specifically for NAR's most compact M-FAK kits
- Mil-Spec hook & loop and 3 heavy duty snaps to ensure kit security

NAR DUAL DIRECTION QUICK DETACH PANEL (4.5 IN. X 6.5 IN.)



DIMENSIONS

- W 4.5 in. x H 6.5 in. x D 0.5 in.
- Weight: 4.5 oz

ITEM #	COLOR
80-0895	Black
80-0896	Coyote

FEATURES/SPECIFICATIONS

- NAR's medium size detach panel with dual direction mount measuring (W 4.5 x H 6.5 in.)
- Utilizes 3 rows by 3 columns vertical mount or 2 rows by 4 columns horizontal mount
- Designed for compatibility with all NAR IFAKs such as horizontal mount for M-FAK and vertical or horizontal mount of the TORK
- Mil-Spec hook & loop and 3 heavy duty snaps to ensure kit security



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277