

K9 BARK

BATTLEFIELD AID & RECOVERY KIT

Mission Ready. Engineered to Deliver.

The K9 BARK™ (Battlefield Aid & Recovery Kit) is a mission-ready canine trauma system engineered to deliver comprehensive point-of-injury care in accordance with the Joint Trauma System (JTS) K9 Tactical Combat Casualty Care (K9TCCC) Guidelines. Purpose built as a hybrid of the proven NAR Mini Medic and NAR Squad Kit Gen 2 platforms, the K9 BARK combines a compact, streamlined design with expanded treatment capability to meet the demands of military, law enforcement, and tactical working dog operations.

Designed for intuitive access under stress, the upper compartment is dedicated to advanced Tier 3/4 interventions, while the lower compartment contains critical Tier 1/2 equipment for immediate lifesaving care. This deliberate organization enables providers to rapidly identify and deploy the appropriate equipment throughout every phase of K9 casualty management, improving efficiency when seconds matter.

From initial lifesaving interventions to advanced battlefield treatment and casualty recovery, the K9 BARK provides operators, handlers, and medics with a single, scalable solution for tactical K9 medical care. Whether supporting military working dogs, law enforcement K9 teams, or search and rescue canines, the K9 BARK delivers the capability, organization, and confidence needed to care for America's four-legged partners in the most demanding operational environments.

FEATURES/SPECIFICATIONS

- Compact, multi-purpose trauma aid bag designed to meet the demands of military, law enforcement, and tactical working dog operations
- Purpose built as a hybrid of the proven NAR Mini Medic and NAR Squad Kit Gen 2 platforms
- Upper compartment is dedicated to advanced Tier 3/4 interventions, while the lower compartment contains critical Tier 1/2 equipment for immediate lifesaving care
- Mini Medic Bag external storage options include 2 pockets and extra drop-down pouch
- MOLLE capability on the sides, bottom and inside back of each bag
- Internal pouch storage of Mini Medic Bag includes 3 clear tear-away, 2 main compartment mesh pouches and 1 larger inner lid pouch
- Squad Gen 2 Bag side straps can be adjusted to contour to body
- Squad Gen 2 Bag can carry front fanny-style or backside
- Durable 500D Nylon & Laser Cut Laminate Molle
- Includes heavy-duty buckle hardware and self-repairing zippers with silent 550 cord pulls
- Sand & water resistant

Learn more about this and other Products with a Mission® at www.NARescue.com



K9 BARK

ITEM #	COLOR	NSN
80-1999	Black	----
80-2000	Coyote	----
80-2001	RGR Green	----

* Explore more K9 Kits at www.NARescue.com

DIMENSIONS

- Mini Medic Bag Dimensions: H 14 in. x W 11 in. x D 5 in. (Depth expandable to 7.25 in.)
- Squad Bag Gen 2 Dimensions: H 7 in. x W 10.5 in. x D 8 in. (expanded D 10 in.)
- Total Weight: 14.9 lbs.



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

K9 BARK

BATTLEFIELD AID & RECOVERY KIT

all KITS INCLUDE:

 QTY 4 Combat Gauze® Z-Fold	 QTY 4 NAR® Gauze, Wound Packing	 QTY 4 ETD® 4 in. Dressing, Flat	 QTY 1 Tactical Suction Device	 QTY 1 BAC-Pack Bougie-Aided Cric	 QTY 1 Handle, Laryngoscope	 QTY 1 #3 Blade, Miller
 QTY 1 Pocket Exam Light	 QTY 2 Alkaline AAA Batteries	 QTY 1 Nylon Pouch	 QTY 2 HyFin® Vent Compact, Twin Pk	 QTY 4 Petrolatum Gauze	 QTY 2 Enhanced ARS Needle, 14G	 QTY 1 Cyclone® Pocket BVM
 QTY 2 Needleless Saline Lock Kit	 QTY 2 Baxter IV Tubing	 QTY 1 NAR Ready Heat 4-Cell	 QTY 1 Emergency Survival Wrap	 QTY 4 Triangular Bandage - ODG	 QTY 2 Dry Sterile Burn Dressing	 QTY 2 SAM® Splint II
 QTY 4 6" Elastic Wrap Bandage	 QTY 2 Surgical Tape, 2 in.	 QTY 2 Surgical Tape, 1 in.	 QTY 1 4" Elastic Adhesive Tape	 QTY 2 1" Coban Wrap Bandage	 QTY 2 Adhesive Tape Bandage	 QTY 1 Tactical Muzzle
 QTY 1 NAR® K-9 Litter	 QTY 1 Dual Head Stethoscope	 QTY 1 Splinter Forceps, 4.5"	 QTY 1 Rectal Thermometer	 QTY 2 Prep Razor	 QTY 1 Sharps Shuttle	 QTY 1 NAR® Mini Duct Tape, 2 in.
 QTY 6 Writeable Name Tape	 PAIR 4 Nitrile Gloves, Lg	 QTY 2 NAR® Trauma Shears, Lg	 QTY 4 K-9 Combat Casualty Card	 QTY 2 Permanent Marker, Large		



K9 BARK

ITEM #	COLOR	NSN
80-1999	Black	----
80-2000	Coyote	----
80-2001	RGR Green	----

* Explore more K9 Kits at www.NARescue.com

Learn more about this and other Products with a Mission® at www.NARescue.com



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277