

HEINE 3000 OPHTHALMOSCOPE FOR CORPSMAN KIT

The Heine® 3000 Ophthalmoscope is designed for examining the interior structures of the eye, especially the retina, in a mini, pocket-size format. This critical high quality optic tool offers all of the important features of a high-end ophthalmoscope combined in a compact, durable, portable configuration.



FEATURES/SPECIFICATIONS

- LED concentrated white light for perfect illumination & crisp retinal detail
- -20 to +20 with 18 lenses; +/-1 to -4/+4 in single steps
- 2-piece head & handle with chrome finish upper section/ refined plastic; Shockproof, sturdy, non-slip
- 5 apertures for full diagnostic capability
- Easy to maintain, flexible
- Integral pocket clip with auto off
- Takes 2 AA batteries (included)

DIMENSIONS

- H 6.75 in. x W 1.25 in. x D 0.9 in.
- Weight: 3.7 oz

HEINE 3000 OPHTHALMOSCOPE

ITEM #	NSN
20-0007	6540-01-491-6651

* Find more medical examination tools at NARescue.com

Learn more about this and other *Products with a Mission®* at www.NARescue.com



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277