

NASOPHARYNGEAL AIRWAY WITH LUBRICANT

The Nasopharyngeal Airway 28F with Lubricant shows North American Rescue's attention to details. This NPA was designed to be dependable, efficient, comfortable, and secure.

Efficient. The product includes a conveniently attached sealed packet of surgical lubricant to facilitate easy insertion.

Designed to be more compatible with human anatomy than the standard oropharyngeal airway, this product is recommended for use as a nasal airway adjunct in semi-conscious or unconscious patients with an intact gag reflex.



Learn more about this and other Products with a Mission® at www.NARescue.com



Made in America

FEATURES/SPECIFICATIONS

- Designed for semi-conscious or unconscious patients with an intact gag reflex
- Packaged with surgical lubricant attached

DIMENSIONS

- Packaged: L 6 in. x W 16 in. x D 0.5 in.
- Weight: 0.5 oz
- Device Length: 117 mm

NASOPHARYNGEAL AIRWAY W/ LUBE

ITEM #	SIZE	LENGTH	NSN
ZZ-0034	28F	117 mm	6515-01-529-1187



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277