

PABC Low Profile Metal Wall Cabinet

An Ideal Solution To Safely Store & Protect

The PABC Low Profile Metal Wall Cabinet easily mounts on the wall with a slim 4 inch clearance. This compact slimline storage solution provides rapid public access to medical equipment when needed. Its white powder painted corrosive resistant galvanized steel construction features a stainless steel continuous hinge and magnetic closure providing quick easy access with one pull to open - No twist and pull. Cases are water resistant. To best protect your critical life savings products cases are best mounted under cover and out of the direct elements. Also included are highly visible "Bleeding Control Kit" stickers for fast identification and convenient "Instructions for Use" located on the interior door panel.



Learn more about this and other Products with a Mission® at www.NARescue.com



FEATURES/SPECIFICATIONS

- Low profile metal cabinet mounts on the wall with a slim 4 inch clearance
- Provides convenient slimline storage for medical equipment
- White powder painted corrosive resistant galvanized steel cabinet with stainless steel continuous hinge
- Magnetic closure provides quick easy access with one pull to open - No twist and pull
- Includes highly visible "Bleeding Control Kit" lettering for fast identification and convenient "Instructions for Use" on the interior door panel

LOW PROFILE CABINET

ITEM #

ZZ-1199

* Explore Public Access Bleeding Control Stations and other cabinet options at www.NARescue.com

DIMENSIONS

- H 15.5 in. x W 22.5 in. x D 4 in.
- Weight: 9.365 lb



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277