

PABC LOW PROFILE METAL STATION - 8 PACK

Save lives! Stop life-threatening bleeding in the first critical minutes!

The first minutes following a traumatic injury are crucial in saving a life. Uncontrolled bleeding can result in death within minutes—even before emergency medical services or other emergency responders can arrive. North American Rescue® (NAR) Public Access Bleeding Control (PABC) Stations provide bystanders and initial first responders with quick and easy access to essential medical equipment to stop life-threatening bleeding in public buildings and populated public areas. This PABC Station includes a low profile metal wall cabinet equipped with eight (8) Individual Bleeding Control Kits (Nylon or Vacuum Sealed) stored in two (2) removable 4 Pocket Mesh Storage Panels with nylon pull handle. Designed to be mounted on a wall at a designated location or to be co-located with Automatic External Defibrillators (AEDs), NAR Public Access Bleeding Control Stations provide publicly accessible medical supplies to control major bleeding.

The low profile metal wall cabinet with 4 inch clearance provides convenient slimline storage for medical equipment. Its white powder painted corrosive resistant galvanized steel construction features a stainless steel continuous hinge and magnetic closure providing quick easy access with one pull to open - No twist and pull. Also included are highly visible "Bleeding Control Kit" stickers for fast identification and convenient "Instructions for Use" located on the interior door panel.

NAR Low Profile Public Access Bleeding Control (PABC) Metal Stations are available in 4 configurations: Basic, Intermediate, Advanced BCD, and Advanced Celox™.

BASIC:

The *Basic version of the NAR PABC Station* provides 8 Individual Bleeding Control Kits (nylon or vacuum sealed) containing lifesaving equipment such as tourniquets, pressure dressings and gauze bandages made publicly accessible to help control bleeding and save lives.

INTERMEDIATE:

Includes all equipment found in the *Basic PABC Station* with each of the 8 Individual Bleeding Control Kits also containing the added capability of treating penetrating chest wounds with the HyFin® Vent Compact Chest Seal Twin Pack (for entrance and exit wounds).

ADVANCED / ADVANCED BCD / ADVANCED CELOX™:

Includes all equipment found in the *Intermediate PABC Station* with each of the complete ruggedized vacuum sealed or nylon pouches also containing the added capability of a hemostatic dressing to enhance blood clotting abilities for faster bleeding control.

PABC STATIONS LOW PROFILE METAL WALL STATIONS

ITEM #	LEVEL	BAG
80-1083	Basic	Vacuum
80-1084	Intermediate	Vacuum
80-1085	Adv. w. BCD	Vacuum
80-1591	Adv. w. CELOX	Vacuum

ITEM #	LEVEL	BAG
80-1080	Basic	Nylon
80-1081	Intermediate	Nylon
80-1082	Adv. w. BCD	Nylon
80-1592	Adv. w. CELOX	Nylon

Learn more about this and other **Products with a Mission®** at www.NARescue.com



FEATURES/SPECIFICATIONS

- Designed to provide bystanders and initial first responders in public buildings and populated public areas with quick, easy access to essential medical equipment for stopping life-threatening bleeding
- Includes eight (8) Individual Bleeding Control Kits packaged in either red nylon bags or clear vacuum sealed pouches and stored in two (2) easily removable Mesh Panels each equipped with 4 pockets and nylon pull handle

- Low profile metal wall cabinet features a slim 4 inch wall clearance with easy pull-to-open magnetic closure for rapid public access to major bleeding control supplies
- Featuring highly visible "Bleeding Control Kit" lettering for fast identification and convenient "Instructions for Use" on the interior door panel
- Stations are available in Basic, Intermediate, Advanced BCD, and Advanced Celox™ care configurations

DIMENSIONS

- H 15.5 in. x W 22.5 in. x D 4 in.
- Station Weight: 18.69 lbs. (*Basic is 17.49 lb*)
- Metal Cabinet Weight: 9.365 lb

* View other Public Access Bleeding Control products at www.NARescue.com



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

PABC LOW PROFILE METAL STATION - 8 PACK



all 8 PACK STATIONS INCLUDE:



all INDIVIDUAL KITS INCLUDE:



+plus, by kit type:



Learn more about this and other **Products with a Mission®** at www.NARescue.com



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277