

PUBLIC ACCESS BLEEDING CONTROL STATIONS

Metal Semi-Recessed: 6 pack & 8 pack

Save lives! Stop life-threatening bleeding in the first critical minutes!

The first minutes following a traumatic injury are crucial in saving a life. Uncontrolled bleeding can result in death within minutes—even before emergency medical services or other emergency responders can arrive. North American Rescue® (NAR) Public Access Bleeding Control (PABC) Stations provide bystanders and initial first responders with quick and easy access to essential medical equipment to stop life-threatening bleeding in public buildings and populated public areas. This PABC Station includes a metal, semi-recessed case equipped with either a six pack (6 vacuum sealed) or eight pack (8 nylon or vacuum sealed) of Individual Bleeding Control Kits contained in an easy-to-open nylon carrying case. Also included is one (1) NAR Responder QuikLitter™ for the 8 pack, or one (1) NAR Quiklitter™ LITE for the 6 pack, for moving injured casualties. Designed to be mounted on a wall at a designated location or to be co-located with Automatic External Defibrillators (AEDs), NAR Public Access Bleeding Control Stations provide publicly accessible medical supplies to control major bleeding.

BASIC:

The *Basic version of the NAR PABC Station* provides 6 or 8 complete kits containing life saving bleeding control equipment such as tourniquets, pressure dressings and gauze bandages to put into the hands of the public to help control bleeding and save lives. It also includes a litter for moving injured casualties.

INTERMEDIATE:

Includes all equipment found in the *Basic PABC Station* with each of the complete kits also containing the added capability of treating penetrating chest wounds with the HyFin® Vent Compact Chest Seal Twin Pack (for entrance and exit wounds).

ADVANCED / ADVANCED BCD / ADVANCED CELOX™:

Includes all equipment found in the *Intermediate PABC Station* with each of the complete ruggedized vacuum sealed or nylon pouches also containing the added capability of a hemostatic dressing to enhance blood clotting abilities for faster bleeding control.

FEATURES/SPECIFICATIONS

- Designed to provide bystanders and initial first responders in public buildings and populated public areas with quick, easy access to essential medical equipment for stopping life-threatening bleeding.
- Each Basic Individual Bleeding Control Kit contains intuitive and easy-to-use tools proven to help to save lives such as the most effective, easiest to use, fastest to apply C-A-T® Tourniquet
- Includes easy-to-open nylon carrying case containing 6 or 8 Individual Bleeding Control Kits packaged in ruggedized, vacuum sealed or nylon pouches and one (1) NAR QuikLitter™ for moving injured casualties
- Available in Basic, Intermediate, Advanced, Adv. BCD, and Adv. Celox™
- Alarm available upon special request. Contact Customer Service for details

Learn more about this and other **Products with a Mission®** at www.NARescue.com



PABC STATIONS - SEMI-RECESSED 8-PACK

ITEM #	LEVEL	#KITS	BAG	WGT
80-0861	Basic	8	Vacuum	19.8 lb
80-0862	Intermediate	8	Vacuum	20.6 lb
80-0864	Adv.	8	Vacuum	21 lb
80-0863	Adv. w. BCD	8	Vacuum	21 lb
80-1578	Adv. w. CELOX	8	Vacuum	21 lb
80-0865	Basic	8	Nylon	21.5 lb
80-0866	Intermediate	8	Nylon	22.3 lb
80-0868	Adv.	8	Nylon	22.7 lb
80-0867	Adv. w. BCD	8	Nylon	22.7 lb
80-1579	Adv. w. CELOX	8	Nylon	22.7 lb

6-PACK

ITEM #	LEVEL	#KITS	BAG	WGT
80-0577	Basic	8	Vacuum	18.7 lb
80-0578	Intermediate	8	Vacuum	19.5 lb
80-0579	Adv.	8	Vacuum	19.5 lb
80-0580	Adv. w. BCD	8	Vacuum	19.5 lb
80-1577	Adv. w. CEL	8	Vacuum	19.5 lb

SEMI-RECESSED METAL CASE FEATURES

- Mounts permanently inside the wall with a low profile, 3.25 in. trim
- Aluminum cabinet with a white finish and clear plexiglass window
- Concealed hinges, magnetic closure & recessed handle designed to open with one quick pull

DIMENSIONS

- Exterior: H 16 in. x L 17 in. x D 3.25 in.
- Interior: H 14 in. x W 15 in. x D 7.25 in.
- Weight: *see table*



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

* View other Public Access Bleeding Control products at www.NARescue.com

PUBLIC ACCESS BLEEDING CONTROL STATIONS

Metal Semi-Recessed: 6 pack & 8 pack



Learn more about this and other **Products with a Mission®** at www.NARescue.com



all 8 PACK STATIONS INCLUDE:



all 6 PACK STATIONS INCLUDE:



all INDIVIDUAL KITS INCLUDE:



+plus, by kit type:



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277