

PUBLIC ACCESS BLEEDING CONTROL

6-PACK - VACUUM SEALED

Save lives! Stop life-threatening bleeding in the first critical minutes!

The first minutes following a traumatic injury are crucial in saving a life. Uncontrolled bleeding can result in death within minutes—even before emergency medical services or other emergency responders can arrive. The North American Rescue® (NAR) Public Access Bleeding Control (PABC) 6-Pack is designed to provide bystanders and initial first responders with quick and easy access to essential medical equipment for stopping life-threatening bleeding. The 6-Pack includes six (6) Individual Bleeding Control Kits that are contained inside of an easy-to-open carrying case. Also included is one (1) NAR Responder QuikLitter™ used for moving injured casualties.

Each of the six Individual Bleeding Control Kits are packaged in a ruggedized vacuum sealed pouch and contain the most effective, easiest to use and fastest to apply C-A-T® Tourniquet, a Responder Emergency Trauma Dressing to be used as a pressure bandage, NAR Wound Packing Gauze for bandaging wounds and step-by-step illustrated pictorial instructions.

The NAR PABC 6-Packs are available in Basic, Intermediate, Advanced, Advanced BCD, and Advanced Celox™ configurations.

BASIC:

The Basic version of the NAR PABC 6-Pack provides 6 complete kits containing life saving bleeding control equipment such as tourniquets, pressure dressings and gauze bandages and puts them into the hands of the public to help control bleeding and save lives. It also includes the NAR Responder QuikLitter™ for moving injured casualties.

INTERMEDIATE:

Includes all equipment found in the Basic Public Access Bleeding Control Pack with each of the 6 complete kits also containing the added capability to treat penetrating chest wounds with the HyFin® Vent Compact Chest Seal Twin Pack (for entrance and exit wounds).

ADVANCED / ADVANCED BCD / ADVANCED CELOX™:

Includes all equipment found in the Intermediate Public Access Bleeding Control Pack with each of the 6 complete kits also containing the added capability of a hemostatic dressing to enhance blood clotting abilities for faster bleeding control.

DIMENSIONS

6-Pack:

- H 7.5 in. x W 13 in. x D 8 in.
- Weight: 5.7 - 6.5 lb

Individual Bleeding Control Kits:

- H 7 in. x W 5 in. x D 2.5 in.
- Weight: 14-16.4 oz

Learn more about this and other Products with a Mission® at www.NARescue.com



FEATURES/SPECIFICATIONS

- Designed to provide bystanders and initial first responders with quick and easy access to essential medical equipment for stopping life-threatening bleeding
- Includes an easy-to-open nylon zippered tote bag containing six (6) Vacuum Sealed Individual Bleeding Control Kits packaged in ruggedized vacuum sealed pouches and one (1) NAR QuikLitter™ Lite for moving injured casualties
- Each Basic Individual Bleeding Control Kit contains intuitive and easy-to-use tools proven to help to save lives such as the most effective, easiest to use, fastest to apply C-A-T® Tourniquet
- Available in Basic, Intermediate, Advanced, Advanced BCD, and Advanced Celox™ care configurations
- Compliant with California Assembly Bill AB2260

PABC 6PK - VAC

ITEM #	OPTION
80-0573	Basic
80-0574	Intermediate
80-0575	Advanced
80-0576	Advanced with BCD
80-1571	Advanced with CELOX™

* Explore other Public Access Bleeding Control options available at NARescue.com



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

PUBLIC ACCESS BLEEDING CONTROL

6-PACK - VACUUM SEALED

all 6-PACKS INCLUDE:



all INDIVIDUAL KITS INCLUDE:



plus, by kit type:



Learn more about this and other Products with a Mission® at www.NARescue.com



6 PACK CONTENTS

- 1 x NAR® Red Nylon Carrying Case
- 1 x NAR® ELITE Litter™ - Responder *LITE*
- 6 x Individual Bleeding Control Kits (Vac)

KIT CONTENTS

- 1 x C-A-T® Tourniquet
- 1 x 6 in. NAR® Responder ETD (Emergency Trauma Dressing)
- 2pr Responder Nitrile Gloves, Lg
- 1 x Responder Trauma Shears, 7.25 in.
- 1 x NAR® Survival Blanket
- 1 x Permanent Marker, Small
- 1 x Just in Time Instruction Card

plus, by kit type:

- **Basic**
 - » 2 x NAR® Wound Packing Gauze
- **Intermediate**
 - » 2 x NAR® Wound Packing Gauze
 - » 1 x HyFin® Compact, Vent Twin Pack
- **Adv.**
 - » 1 x Combat Gauze LE, Z-Fold Hemostatic Bandage
 - » 1 x NAR® Wound Packing Gauze
 - » 1 x HyFin® Compact, Vent Twin Pack
- **Adv. BCD**
 - » 1 x QuikClot Bleeding Control Dressing
 - » 1 x NAR® Wound Packing Gauze
 - » 1 x HyFin® Compact, Vent Twin Pack
- **Adv. CELOX**
 - » 1 x 5 ft. Celox™ First Aid Z-Fold Hemostatic Dressing
 - » 1 x NAR® Wound Packing Gauze
 - » 1 x HyFin® Compact, Vent Twin Pack

PABC 6PK - VAC

ITEM #	OPTION
80-0573	Basic
80-0574	Intermediate
80-0575	Advanced
80-0576	Advanced with BCD
80-1571	Advanced with CELOX™

* Explore other Public Access Bleeding Control options available at NARescue.com



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277