

# PUBLIC ACCESS BLEEDING CONTROL STATIONS

Clear Station: 6 pack & 8 pack

## Save lives! Stop life-threatening bleeding in the first critical minutes!

The first minutes following a traumatic injury are crucial in saving a life. Uncontrolled bleeding can result in death within minutes—even before emergency medical services or other emergency responders can arrive. North American Rescue® (NAR) Public Access Bleeding Control (PABC) Stations provide bystanders and initial first responders with quick and easy access to essential medical equipment to stop life-threatening bleeding in public buildings and populated public areas. This PABC Station includes a clear wall mounted storage case equipped with six (6) Basic Vacuum Sealed Individual Bleeding Control Kits contained in an easy-to-open nylon zippered tote bag. Also included is one (1) NAR® ELITE Litter™ - Responder LITE for moving injured casualties. Designed to be mounted on a wall at a designated location or to be co-located with Automatic External Defibrillators (AEDs), the NAR Public Access Bleeding Control Station provides publicly accessible medical supplies to control major bleeding.

### BASIC:

The Basic version of the NAR PABC Station provides 6 complete ruggedized vacuum sealed pouches containing life saving bleeding control equipment such as tourniquets, pressure dressings and gauze bandages to put into the hands of the public to help control bleeding and save lives. It also includes the NAR® ELITE Litter™ - Responder LITE for moving injured casualties.

### INTERMEDIATE:

Includes all equipment found in the Basic PABC Station with each of the 6 complete ruggedized vacuum sealed pouches also containing the added capability of treating penetrating chest wounds with the HyFin® Vent Compact Chest Seal Twin Pack (for entrance and exit wounds).

### ADVANCED / ADVANCED BCD / ADVANCED CELOX™:

Includes all equipment found in the Intermediate PABC Station with each of the 6 complete ruggedized vacuum sealed pouches also containing the added capability of a hemostatic dressing to enhance blood clotting abilities for faster bleeding control.

## FEATURES/SPECIFICATIONS

- Designed to provide bystanders and initial first responders in public buildings and populated public areas with quick, easy access to essential medical equipment for stopping life-threatening bleeding.
- Includes easy-to-open nylon zippered tote bag containing six (6) Individual Bleeding Control Kits packaged in ruggedized vacuum sealed pouches and one (1) NAR® ELITE Litter™ - Responder LITE for moving injured casualties
- Alarm available upon special request. Contact Customer Service for details
- Each Basic Individual Bleeding Control Kit contains intuitive and easy-to-use tools proven to help to save lives such as the most effective, easiest to use, fastest to apply C-A-T® Tourniquet
- The rugged wall mounted clear storage case provides high visibility access as well as physical and environmental protection of medical supplies
- Stations are available in Basic, Intermediate, Advanced, Advanced BCD, and Advanced Celox™ care configurations
- Compliant with California Assembly Bill AB2260

Learn more about this and other **Products with a Mission®** at [www.NARescue.com](http://www.NARescue.com)



## PABC STATIONS - CLEAR 8-PACK

| ITEM #  | LEVEL         | #KITS | BAG    | WGT      |
|---------|---------------|-------|--------|----------|
| 80-0471 | Basic         | 8     | Vacuum | 21.8 lb  |
| 80-0472 | Intermediate  | 8     | Vacuum | 22.6 lb  |
| 80-0473 | Adv.          | 8     | Vacuum | 23 lb    |
| 80-0526 | Adv. w. BCD   | 8     | Vacuum | 23 lb    |
| 80-1575 | Adv. w. CELOX | 8     | Vacuum | 23 lb    |
| 80-0474 | Basic         | 8     | Nylon  | 23.55 lb |
| 80-0475 | Intermediate  | 8     | Nylon  | 24.3 lb  |
| 80-0476 | Adv.          | 8     | Nylon  | 24.7 lb  |
| 80-0527 | Adv. w. BCD   | 8     | Nylon  | 23 lb    |
| 80-1576 | Adv. w. CELOX | 8     | Nylon  | 23 lb    |

## 6-PACK

| ITEM #  | LEVEL        | #KITS | BAG    | WGT     |
|---------|--------------|-------|--------|---------|
| 80-0581 | Basic        | 6     | Vacuum | 20.7 lb |
| 80-0582 | Intermediate | 6     | Vacuum | 21.5 lb |
| 80-0583 | Adv.         | 6     | Vacuum | 21.5 lb |
| 80-0584 | Adv. w. BCD  | 6     | Vacuum | 21.5 lb |
| 80-1574 | Adv. w. CEL  | 6     | Vacuum | 19.5 lb |

## CLEAR POLYCARBONATE CASE FEATURES

- Visibility at any angle
- External mounting tabs
- Includes a wire shelf that supports the Public Access Bleeding Control Pack
- Exterior has rounded corners to help prevent injury
- Comes complete with the mounting hardware

## DIMENSIONS

- Exterior: H 18 in. x W 15.5 in. x D 10 in.
- Weight: *see table*

\* View other Public Access Bleeding Control products at [www.NARescue.com](http://www.NARescue.com)



**NORTH AMERICAN RESCUE®**  
[www.NARescue.com](http://www.NARescue.com) • 888.689.6277

# PUBLIC ACCESS BLEEDING CONTROL STATIONS

*Clear Station: 6 pack & 8 pack*



6 Pack



6 Pack  
Vacuum Sealed



6 Pack Nylon

Learn more about this and  
other **Products with a Mission®**  
at [www.NARescue.com](http://www.NARescue.com)



## all 8 PACK STATIONS INCLUDE:




## all 6 PACK STATIONS INCLUDE:



## all INDIVIDUAL KITS INCLUDE:



+plus, by kit type:



**NORTH AMERICAN RESCUE®**  
[www.NARescue.com](http://www.NARescue.com) • 888.689.6277