


HYFIN® CHEST SEAL DRESSING - COMBO PACK

The HyFin® Chest Seal Combo Pack from North American Rescue provides one HyFin® Vented Chest Seal and one HyFin® Occlusive Chest Seal packaged in a more compact, rugged, easy-to-open twin foil pouch Combo Pack. Each HyFin® chest seal is individually packaged within the Combo Kit and includes a gauze pad to wipe the wound surface prior to application. The two different seals are each individually packaged and are designed to allow the user to apply a seal to both an entry / exit wound while also giving the option to only apply one as needed and store the other as the situation requires.

Watch the HyFin®
Occlusive Chest Seal
Overview Video



Watch the HyFin® Vent
Chest Seal Illustrated
Instructions



Learn more about this and
other Products with a Mission®
at www.NARescue.com



Includes one HyFin®
Vented Chest Seal &
one HyFin® Occlusive
Chest Seal

FEATURES/SPECIFICATIONS

- Versatile Combo Pack designed for the treatment of both entry/exit or multiple penetrating injuries to the chest
- Featuring NAR's patented HyFin® Chest Seals with enhanced adhesive technology designed for a superior seal in the most adverse conditions, including sweaty or hairy casualties
- Each HyFin® measures 6-inch x 6-inch and contains a gauze pad to wipe the wound surface prior to application
- Includes one HyFin® Vented Chest Seal & one HyFin® Occlusive Chest Seal individually packaged in clearly marked and easy-to-open twin foil pouch
- Large, easily-visible red elliptical tabs designed to allow for single step peel-and-apply application and the burping of a wound if necessary

HYFIN® COMBO PACK

ITEM #	NSN
10-0073	6510-01-654-6787

* Explore more HyFin® chest seal options at www.NARescue.com

DIMENSIONS

- Device: H 6 in. x W 6 in. (excluding tab)
- Packaged Folded: H 7.5 in. x W 4.5 in. x D 0.25 in.
- Packaged Unfolded: L 7.5 in. x W 9 in. x D 0.13 in.
- Weight: 2.5 oz



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277