

CELOX™ Z-FOLD HEMOSTATIC GAUZE

5FT. CELOX Z-FOLD HEMOSTATIC GAUZE

Learn more about this and other Products with a Mission® at www.NARescue.com



CELOX™ Gauze in a versatile Z-fold configuration

Celox™ Gauze works independently of the body's clotting mechanism. Users apply the gauze directly on to the bleeding source and hold pressure for three minutes to stop bleeding from arterial injuries, gunshot wounds, road traffic accidents and other moderate to severe bleeding injuries.

Celox™ blood clotting products have very high surface area chitosan flakes. When they come in contact with blood, Celox™ swells, gels, and sticks together to make a gel like plug. Celox™ blood clotting agent does not set off the normal clotting cascade, it only clots the blood it comes directly into contact with.

Celox™ Gauze hemostat has been CoTCCC approved for DoW use.



FEATURES/SPECIFICATIONS

- Stops life-threatening bleeding found in severe traumatic injuries
- Works independently of the body's normal clotting mechanisms

CELOX™ Z-FOLD

ITEM #	NSN
30-0064	6510-01-623-9910

* Explore more hemostatic products at www.NARescue.com



DIMENSIONS

- Packaged: H 5.25 in. x W 4.25 in. x D 0.5 in.
- Gauze: W 3 in. x L 5 ft
- Weight: 11 oz



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277