

QINFLOW WARRIOR LITE INTRO SET

LITE or EXTRA POWER BATTERY

High-performing yet compact and lightweight, the Warrior Lite is ideally configured for rescue teams operating under space and weight constraints, as well as first responders, emergency medical services (EMS), helicopter emergency medical services (HEMS), and critical care transport teams.

The Warrior Lite offers two battery options to accommodate different operational needs. The large ("Extra Power") battery is designed to warm up to 6 units of cold-stored whole blood or up to 4.5 liters of room-temperature fluids. The small ("Lite") battery provides approximately half that capacity. Both battery options support rapid warming, delivering a cold-stored whole blood unit in less than 3 minutes.

The Extra Power Battery is well-suited for scenarios requiring higher volumes of blood or fluid warming, such as prolonged transports, consecutive missions, and incidents involving multiple critically injured patients, including emergency department settings. The Lite battery is the preferred choice when minimizing size and weight is the primary consideration.



FEATURES/SPECIFICATIONS

- High-performing, yet compact and lightweight
- Ideally configured for rescue teams operating under space and weight constraints
- Supports rapid warming, delivering a cold-stored whole blood unit in less than 3 minutes
- "Lite" Battery option is the preferred choice when minimizing size and weight is the primary consideration

QINFLOW WARRIOR LITE INTRO SET

ITEM #	OPTION	NSN
35-0046	Lite Battery	----
35-0071	Extra Power Battery	----

R_X ONLY

Learn more about this and other Products with a Mission® at www.NARescue.com



DIMENSIONS

Lite Battery:

- Dimensions: H 3.3 in. x W 3.46 in. x L 4.52 in.
- Weight: 1.76 lbs.

Extra Power Battery:

- Dimensions: H 4.13 in. x W 3.46 in. x L 4.52 in.
- Weight: 2.39 lbs.

* Explore more QinFlow Products at www.narescue.com



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

QINFLOW WARRIOR LITE INTRO SET

LITE or EXTRA POWER BATTERY

LITE KITS INCLUDE:



EXTRA POWER KITS INCLUDE:



**Compact
Disposable Unit
Sold Separately**



KIT CONTENTS

- 1 x QinFlow Warrior Lite
- 1 x Li-Ion Lite - OR - Extra Power Battery
- 1 x Lite Battery Charger - OR - Extra Power Battery

Learn more about this and
other Products with a Mission®
at www.NARescue.com



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277