

Learn more about this and other Products with a Mission® at www.NARescue.com



QINFLOW WARRIOR LITE INTRO SET

LITE or EXTRA POWER BATTERY

Box of 12 consumable single-patient sterile Compact Disposable Units (CDUs). The CDU is the warming component with up to 3 years of shelf life and it is compatible with all Warrior and Warrior lite models.

NOTE: The CDU must be used together with one of the **Warrior Lite Intro Sets** (Lite Battery or Extra Power Battery).

The CDU can be easily transferred from one Warrior model to the next to facilitate a complete continuum of care implementation; thus, simplifying patient handoffs and reducing the overall cost of ownership by requiring fewer consumables.



QinFlow Warrior Lite Intro Set Sold Separately



FEATURES/SPECIFICATIONS

- Sterile, consumable, single-patient use
- Comes in a box of 12
- Warming component with up to 3 years of shelf life
- Compatible with all Warrior and Warrior Lite Models

QINFLOW CDU ^{Rx ONLY}

ITEM #	NSN
35-0047	----

* Explore more QinFlow Products at www.narescue.com

DIMENSIONS

- Dimensions: H 2.84 in. x W 4.63 in. x D 2.70 in.
- Weight: 0.26 lb.



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277