

TALON II MODEL 90C

COLLAPSIBLE HANDLE LITTER

The upgraded Talon II 90C is the most versatile and compact evacuation and treatment litter available today! Developed to meet the U.S. Army's requirement to provide urgent casualty evacuation from the battlefield. The new design allows patient transport in restricted compartments such as the up-armored M1114 Humvee and other non-conventional commercial vehicles. By simply extending the ergonomically designed handles it transitions to a standard NATO compatible evacuation platform. The Talon 90C is the ultimate evacuation platform.



FEATURES/SPECIFICATIONS

- Advanced, innovative combat evacuation system featuring patented portable quad-fold design with retractable handles
- Developed to U.S. Army requirements
- Certified Safe-to-Fly U.S. Army Rotary Wing
- Retractable, ergonomically designed handles
- Automatic locking hinges
- Slip-resistant, chemical-resistant ripstop fabric
- 6 IV-attachment points
- 2 patient-securing straps
- U.S. Patents 6,842,923 & 5,598,592, Patent(s) Pending



Learn more about this and other Products with a Mission® at www.NARescue.com



By extending the ergonomically designed handles it transitions to a standard NATO compatible evacuation platform

TALON II 90C

ITEM #	NSN
60-0002	6530-01-504-9051
CE 60-0002-EU	6530-01-504-9051

* Explore other evacuation products at www.NARescue.com

DIMENSIONS

Folded with handles collapsed:

- H 20.5 in. x W 8.5 in. x D 9.5 in.

Open with handles collapsed:

- L 78 in. x W 22.5 in. x D 6 in.

Open with handles extended:

- L 90 in. x W 22.5 in. x D 6 in.

Weight:

- 16 lb 6 oz (7.4 kg)

Static Payload:

- Load Tested to 450lbs (204kg)
- Load Capacity up to 1,500lbs (680kg)

Ground Clearance:

- 15 in.



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277