

RAVEN MODEL 90C

BIFOLD LITTER

Official Litter for the U.S. NAVY.

The challenge for cross-capability between Homeland Security and the Department of Defense was "mission possible" for North American Rescue. We designed the Raven® Model 90C as a cost-effective aid for rescuers faced with the challenges of natural disasters and domestic terrorism.

Our unique solution includes such important features as bi-fold design, and durable construction for reliable long-term storage. The stretcher allows immediate evacuation on DOD aircraft and integration into DOD medical treatment facilities. Importantly, it also works as a treatment table at remote locations in all environments. The Raven® 90C has been tested for chemical warfare agents.

Raven® 90C also has a family of support items such as IV poles, stretcher stands, and backrest.

FEATURES/SPECIFICATIONS

- Certified Safe-to-Fly U.S. Air Force
- Bifold design for reliable long-term storage.
- Ergonomically designed handles
- 2 IV Attachment Points
- 2 Patient Securing Straps
- Slip-resistant, chemical-resistant, rip-stop fabric tested & approved by Federal & private labs using live warfare agents
- U.S. Patents 6,842,923 & 5,598,592, Patent(s) Pending

TALON II 90C

ITEM #	NSN
60-0001	6530-01-432-5114
60-0001-EU	6530-01-432-5114



* Explore other evacuation products at www.NARescue.com



DIMENSIONS

Folded:

- H 45.5 in. x W 7.5 in. x D 8 in.

Open:

- L 90 in. x W 22.5 in. x D 6 in.

Weight:

- 16 lb 8 oz (7.5 kg)

Static Payload:

- Load Tested to 450lbs (204kg)
- Load Capacity up to 1,500lbs (680kg)

Ground Clearance:

- 15 in.



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

Learn more about this and other Products with a Mission® at www.NARescue.com