

T3 LIGHTWEIGHT ASSAULT LITTER

Our Newest Evolution of Evacuation Platforms...Weighing less than 12 lb! The T3 Lightweight Assault Litter is designed to support Special Operation Forces where every pound matters. State of the art materials and technologies have allowed us to reduce the weight without compromising the overall strength and performance of the evacuation platform. The overall length of 81 in. enables floor loading on a UH-60 helicopter while allowing for maximum maneuverability during casualty evacuation in urban terrain. Weighing in at less than twelve pounds, the T3 is the lightest rigid litter in the world.



FEATURES/SPECIFICATIONS

- Weighs less than 12 lb
- Maximum maneuverability during TACEVAC in urban terrains
- Overall length of 81 in. enables floor loading on a UH-60 helicopter
- 2 straps for securing patients
- 4 IV attachment points
- Ergonomically designed handles
- Quad-fold design for maximum cube space storage

Learn more about this and other Products with a Mission® at www.NARescue.com



T3 LITTER

ITEM #	NSN
60-0022	6530-01-663-1222

** Explore other evacuation products at www.NARescue.com*

DIMENSIONS

Folded:

- H 21 in. x W 8.75 in. x D 7.25 in.

Open:

- L 81 in. x W 22.5 in. x D 6 in.

Weight:

- 11 lb 13 oz (5.34 kg)

Static Payload:

- Load Tested to 450lbs (204kg)
- Load Capacity up to 1,500lbs (680kg)

Ground Clearance:

- 1.5 in.



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277