

OSL OVERSIZED LITTER

Designed to safely evacuate larger patients with both comfort and ease. This design has an increased patient surface area, extra padded carrying handles, flame-retardant cover, and additional stirrups offering firmer patient support as well as better weight distribution on the floor. Like the Raven 90C, the OSL has a bi-fold configuration. The cover material is a very durable Nomex/Kevlar blend. Also, a supplemental Back Rest (not shown) to support the patient at both 90° and 45° angles may be purchased separately.

FEATURES/SPECIFICATIONS

- Certified safe-to-fly by the U.S. Air Force
- Flame-retardant cover made of a durable Nomex/Kevlar-blend fabric
- 6 stirrups for firmer patient support and better weight distribution
- Ergonomically designed handles
- Additional carrying loops
- 6 IV attachment points
- 3 patient securing straps
- Attached storage pouch
- Bifold configuration for easy storage
- U.S. Patent 5,598,592, Patent(s) Pending



Learn more about this and other Products with a Mission® at www.NARescue.com



OSL

ITEM #	NSN
60-0023	6530-01-552-7708

* Explore other evacuation products at www.NARescue.com

DIMENSIONS

Folded:

- H 45 in. x W 9 in. x D 8 in.

Open:

- L 90 in. x W 31 in. x D 6 in.

Weight:

- 22 lb 6 oz (10.2 kg)

Static Payload:

- Load Tested to 450lbs (204kg)
- Load Capacity up to 1,500lbs (680kg)

Ground Clearance:

- 15 in.



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277