

# OPIOID OVERDOSE KITS



## Compact and Durable Opioid Antidote Kit...


The Opioid Overdose Kits from North American Rescue are designed to help first responders with the intranasal administration of an opioid antidote to patients with known or suspected opioid overdose in the prehospital setting.

These compact, durable kits are housed in protective carrying cases and include an array of equipment configurations to meet varying local protocols for the emergency treatment of opioid overdose.

Kits are comprised of a pair of nitrile gloves, CPR Microshield, Mucosal Atomizer and come with or without Naloxone Hydrochloride (Narcan®) - 2mg/2ml. Kit components for resupply are also available for individual purchase.

The rugged nylon carrying case features a protective foam insert to protect the medication from being damaged and a built-in eyelet on the outer zipper to facilitate a tamper-resistant security seal for the kit and its contents.

Also included is a front mounted kit identification card to track the expiry and lot information of the medication.



**NORTH AMERICAN RESCUE®**  
[www.NARescue.com](http://www.NARescue.com) • 888.689.6277



## No Medication Opioid Overdose Kit



1 Mucosal Atomizer

1 CPR Face Shield

1 Pair Nitrile Gloves

1 Expiry/Lot MED ID Card

## Single Dose Opioid Overdose Kit



1 Mucosal Atomizer

1 CPR Face Shield

1 Pair Nitrile Gloves

1 Expiry/Lot MED ID Card

1 Narcan<sup>®</sup> Naloxone HCL

## Double Dose Opioid Overdose Kit



2 Mucosal Atomizer

1 CPR Face Shield

2 Pair Nitrile Gloves

1 Expiry/Lot MED ID Card

2 Narcan<sup>®</sup> Naloxone HCL

### Features:

- Compact, durable protective carrying case with equipment to help first responders with the intranasal administration of an opioid antidote to treat opioid overdoses in the prehospital setting
- Rugged nylon carrying case with a protective foam insert to protect the medication from being damaged
- Zipper has a built-in eyelet for the use of a security seal to easily identify if the kit has been opened or tampered with
- Includes a kit identification card to track the expiry and lot information of the medication
- Kits are available with and without Naloxone Hydrochloride (Narcan<sup>®</sup>)

### Dimensions:

- L 6.5 in. x W 3.25 in. x D 2 in.
- Weight: 3.3 oz



*Individual Components and ReSupply Kits are also available, please see [www.NARescue.com](http://www.NARescue.com) for more information.*



**NORTH AMERICAN RESCUE<sup>®</sup>**  
[www.NARescue.com](http://www.NARescue.com) • 888.689.6277

PRODUCTS WITH A  
**MISSION<sup>®</sup>**