

# NAR TRIANGULAR BANDAGE



## **NAR Triangular Bandage (Cravat):**

Part of the Massive Hemorrhage control product line, the NAR Triangular Bandage has a very small footprint - only measuring 3.5 inches at its widest point and weighing only 1.6 ounces.

Coming pre-packaged with two safety pins for ease of use, and available in OD Green, this is a must have for any first aid kit.

## **Features:**

- For the bandaging of wounds, the immobilization of fractures and dislocations, and as a sling for the support of an injured part of the body
- Packaged for low cube space
- Constructed of muslin cloth with two safety pins to aid in bandage configuration and securing
- Easily fits into any first aid kit or pocket

## **Specifications:**

- Packaged: 3.5 in. x 2 in. x .625 in.
- Opened: 37 in. x 37 in. x 52 in.
- Weight: 1.6 oz.



Item#	NSN#
30-0089	6510-00-201-1755



**NORTH AMERICAN RESCUE®**

www.NARescue.com • 888.689.6277

PRODUCTS WITH A  
**MISSION®**